

Avser Artiklar: 45-544, 45-539, 45-541, 45-543



Cannula Cleaning Brushes

These brushes are ideal for cleaning liposuction cannulas, and other general use instruments.

- Designed in a way that could assist within the lumen during the cleaning process
- Made with anti-microbial nylon bristles and stainless steel twisted wire shafts
- Color-coded for easy identification and reordering

Order#	Description	UOM
	Nylon Twisted Wire Brush, 12"/30cm	
45-542	2mm dia, Red Bead	3/pk
45-543	3mm dia, White Bead	3/pk
45-544-2	4mm dia, Aqua Bead	3/pk
45-545-2	5mm dia, Green Bead	3/pk
45-546-2	6mm dia, Black Bead	3/pk
	Nylon Twisted Wire Brush, 16"/40.5cm	
45-539	5mm dia, Green Bead	3/pk
45-540	7mm dia, Yellow Bead	3/pk
45-541	10mm dia, Pink Bead	3/pk



General Instrument Brush Cleaning Set

Order#	Description	UOM
61250	General Instrument Brush Cleaning Set	1ea

Includes:

- | | |
|---|---|
| 1 pk 5ft suction tube cleaning brushes (45-905) | 1 pk Single-ended, fine stainless steel bristle brush (45-303SS) |
| 1 pk 6ft suction tube cleaning brushes (45-906) | 2 pk Single-ended, plastic handle, nylon bristle brush (45-30380) |
| 1 pk 8ft suction tube cleaning brushes (45-906) | 1 pk 3mm dia, nylon bristle brush, 12" (45-543) |
| 1 pk 10ft suction tube cleaning brushes (45-910) | 1 pk 5mm dia, nylon bristle brush, 16" (45-539) |
| 1 pk 12ft suction tube cleaning brushes (45-912) | 1 pk 7mm dia, nylon bristle brush, 16" (45-540) |
| 2 pk Double-ended nylon bristle brush (M-16) | 1 pk 10mm dia, nylon bristle brush, 24" (45-553) |
| 1 pk 9" utility brush, nylon bristle, plastic handle (45-825) | |

Conversion Chart

mm	inch	french
1.0	0.039	3
1.35	0.053	4
1.67	0.066	5
2.0	0.079	6
2.3	0.092	7
2.7	0.105	8
3.0	0.118	9
3.3	0.131	10
3.7	0.144	11
4.0	0.158	12
4.3	0.17	13
4.7	0.184	14
5.0	0.197	15
5.3	0.21	16
5.7	0.223	17
6.0	0.236	19
6.3	0.249	20
6.7	0.263	21
7.3	0.288	22
8.0	0.315	24
8.7	0.341	26
9.3	0.367	28
10.0	0.393	30
10.7	0.419	32
11.3	0.445	34