



# evercare® Surgical hood Anastasia

REF: 2798

evercare® Headweares are designed to help prevent foreign particles such as hair or skin scales falling from medical personnel into the wound during surgical interventions. These very popular surgical hoods come with different special features and designs, such as integrated sweatbands for additional comfort, or open design for easier dressing placed on the head instead of over the head. evercare® Surgical hood Anastasia operating helmet is especially designed for anesthetic staff, glasses users or those with limited hearing capabilities. The thin net over the ears makes it easier to hear each other in a noisy OR, it also makes it easy to use stethoscopes or place glasses. The hood is equipped with an extra-large sweatband for extra comfort. There is shoulder protection, a soft binding round the face and tie-band for individual fitting. Made in a soft and breathable viscose nonwoven material (100% fossile-free.)

## Key Features

- ▶ Reduces the risk of particles falling into the wound
- ▶ Universal size
- ▶ Shoulder protection
- ▶ Integrated large sweatband
- ▶ Specially designed for anesthetic staff
- ▶ Recommended to be used by persons with bad hearing
- ▶ Tie-band




REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2798	Surgical hood Anastasia, shoulder protection, tie-band, large sweatband, Green	Universal	50 / 300

Sales Unit: Transport Carton

# Technical Summary

## Regulations

The product complies with legal requirements for medical devices.



<b>REACH regulation:</b>	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
<b>Biological evaluation:</b>	EN ISO 10993-10
<b>Compliance with legal requirements:</b>	Medical Device Regulation (EU) 2017/745, MDR class I
<b>Quality standards:</b>	EN ISO 13485
<b>Label information and packaging compliance:</b>	EN ISO 15223-1, EN 1041

## Material

Product Part	Raw Material & Characteristics
<b>Main material</b>	Viscose nonwoven, green
<b>Sweatband and face lining</b>	Viscose spunlance nonwoven, white
<b>Elastic band</b>	Nylon and synthetic rubber, latex free, white
<b>Sewing thread</b>	Nylon and spun polyester, white
<b>Tie bands</b>	PP Spunbond, white
<b>Ear cover</b>	PES, white

## Warnings & Recommendations

<b>Warnings</b>	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	3 Years

## Packing information

Package	Material	Size ( L x W x H ) mm
<b>Single Pack</b>	N/A	N/A
<b>Dispenser Pack</b>	Carton, TCF or ECF quality	300 x 185 x 145
<b>Transport Carton</b>	Carton, TCF or ECF quality	385 x 310 x 455
<b>Packing Tape</b>	PVC-free, Transparent	N/A

## Barcodes

REF	Single Pack	Dispenser	Transport Carton
<b>2798</b>	N/A	6438129321630	6438129521634

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.