



evercare® Surgical Hood Jess

REF: 2799

evercare® Headweares are designed to help prevent foreign particles such as hair or skin scales falling from medical personnel into the wound during surgical interventions. These very popular surgical hoods come with different special features and designs, such as integrated sweatbands for additional comfort, or open design for easier dressing placed on the head instead of over the head. evercare® Surgical Hood Jess is designed for users that don't want to pull a helmet over the head. Instead the helmet is placed on the head and closed with a tape strap underneath the chin. The tape makes it easy to adjust the size individually and the user will avoid irritating tie-bands. There is also shoulder protection and an extra-large sweatband to allow for extra comfort. Jess is new in our headwear family and has quickly become popular and is being asked for by medical staff who prefer more flexibility and possibilities for adjustments. Since this hood is without rubber band Jess has also become popular by medical staff who get skin irritation / have sensitive skin. Jess is made in a soft and breathable viscose nonwoven material (100% fossile-free).

Key Features

- ▶ Reduces the risk of particles falling into the wound
- ▶ Universal size
- ▶ Shoulder protection
- ▶ Integrated large sweatband
- ▶ Tape adjustment in the back




REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2799	Surgical hood Jess, shoulder protection, large sweatband, Blue	Universal	50 / 300

Sales Unit: Transport Carton

Technical Summary

Regulations

The product complies with legal requirements for medical devices.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR class I
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN 1041

Material

Product Part	Raw Material & Characteristics
Main material	Viscose nonwoven, blue
Sweatband	Viscose spunlance nonwoven, white
Sewing thread	Nylon and spun polyester, white
Double face tape	Double adhesive tape, latex-free, white

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	3 Years

Packing information

Package	Material	Size (L x W x H) mm
Single Pack	N/A	N/A
Dispenser Pack	Carton, TCF or ECF quality	300 x 185 x 145
Transport Carton	Carton, TCF or ECF quality	385 x 310 x 455
Packing Tape	PVC-free, Transparent	N/A

Barcodes

REF	Single Pack	Dispenser	Transport Carton
2799	N/A	6438129321647	6438129521641

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.