



Product description

Order number	95.1717.011
Product description	Transport case for diagnostic samples in accordance with packaging regulation P 650, type B13, Flat design (LxWxH) 395 x 395 x 300 mm, with two additional loops for sealing, plastic bag with closure clip as secondary packaging, including a

Product characteristics

Design	Flat design for block rack
Example of use	External packaging according to P 650/ ADR
Shape of inner tube	Square

Size

Width of product	395 mm
Length of product	395 mm

Material & colours

Product material	High Density Polyethylen (HD-PE)
------------------	----------------------------------

Purity & certification

Satisfies the requirement	P 650, ADR
Product in accordance with MDD or IVDD	No marking necessary
CE certified	No CE marking
Lot or serial number	non batched

Packaging

Pcs. / inner box	1
Pcs. / box	1
Pcs. / pallet	30
Minimum order qty.	1
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0,0538 cbm
Weight of product	2800 g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification
Transport case B13**

Page 2

Weight of case	3,3 KG
EAN of inner box	4038917133197
EAN of outer box	4038917069250

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2016-05-31

This document was prepared by EDP support and is valid without signature.



SARSTEDT