



Product description

Order number	95.1717.003
Product description	Transport case for diagnostic samples in accordance with packaging regulation P 650, type B 17, high design (length 395mm x width 395mm x height 435mm) with enlarged inner part, plastic bag with closure clip as secondary packaging, incl. 1 superabsor

Product characteristics

Design	Without lateral rim
Example of use	External packaging according to P 650/ ADR
Shape of inner tube	Square

Size

Width of product	395 mm
Length of product	395 mm

Material & colours

Product material	High Density Polyethylen (HD-PE)
------------------	----------------------------------

Purity & certification

Satisfies the requirement	P 650, ADR
Product in accordance with MDD or IVDD	No marking necessary
CE certified	No CE marking
Lot or serial number	non batched

Packaging

Pcs. / inner box	1
Pcs. / box	1
Pcs. / pallet	18
Minimum order qty.	1
Depth of case	403 mm
Width of case	403 mm
Height of case	460 mm
Case volume	0,0747 cbm
Weight of product	4500 g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification
Transport case B19**

Page 2

Weight of case	4,5 KG
EAN of inner box	4038917133173
EAN of outer box	4038917069861

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2016-05-31

This document was prepared by EDP support and is valid without signature.



SARSTEDT