



Product description

Order number	72.383
Product description	CryoPure 5.0 ml cryo tube with white QuickSeal screw cap, with external thread, 4.5 ml sample volume, sterile, non-pyrogenic, non-mutagenic and free from cytotoxins, 25 pieces in the bag, 250 pieces packed in the box and 1.000 pieces packed in the m

Product characteristics

Label/ Print	With printed label
Graduation	Yes
Cap	Screw cap
Internally or externally threaded cap	External thread
Suitable for	93.875.xxx

Size

Sample volume	4,5 ml
Length of product	90,5 mm

Material & colours

Product material	Polypropylene (PP)
Cap material	High Density Polyethylene (HD-PE)
Colour of cap	White

Purity & certification

Product in accordance with MDD or IVDD	In-Vitro Product
CE certified	CE
Sterilisation	Electron irradiation
Pyrogen-free	Pyrogen-free
Lot or serial number	always batched

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Pcs. / inner box	250
Pcs. / box	1000
Pcs. / pallet	52000
Minimum order qty.	250
Depth of case	598
Width of case	378
Height of case	135
Case volume	0,0305 cbm
Weight of case	5,05 KG
EAN of inner box	4038917081412
EAN of outer box	4038917080651

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

