



Product description

Order number	93.853
Product description	Block rack D13 with 50 grooves for Monovettes and tubes up to a diameter of 13 mm, without handles, Grey

Product characteristics

Label/ Print	With printed label
Design	Without handle
Number of wells	50
Storage capacity	50
Autoclavable (max. °C)	121 °C

Size

Diameter of opening	13 mm
Width of product	96 mm
Length of product	180 mm

Material & colours

Product material	Polypropylene (PP)
------------------	--------------------

Purity & certification

Product in accordance with MDD or IVDD	No marking necessary
CE certified	No CE marking
Lot or serial number	not always batched

Packaging

Pcs. / inner box	1
Pcs. / box	16
Pcs. / pallet	1440
Minimum order qty.	1
Depth of case	398
Width of case	368
Height of case	120
Case volume	0,0176 cbm
Weight of case	3,25 KG
EAN of inner box	4038917084390

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





EAN of outer box

4038917074018

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2016-05-31

This document was prepared by EDP support and is valid without signature.



SARSTEDT