



Product description

Order number	72.706.400
Product description	SafeSeal Micro tube 1.5 ml with attached lid with a wide hinge and retaining cams, PP, neutral, PCR Performance Tested certified, 250 pieces packed in the bag

Product characteristics

Graduation	Yes
Cap	SafeSeal
Recommended centrifugation	30000 x g

Size

Sample volume	2 ml
Diameter	10,8 mm
Length of product	39 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	Transparent
Cap material	Polypropylene (PP)
Colour of cap	Transparent

Purity & certification

Product in accordance with MDD or IVDD	No marking necessary
CE certified	No CE marking
Purity standard	PCR Performance Tested Quality
Lot or serial number	always batched

Packaging

Pcs. / inner box	1000
Pcs. / box	2000
Pcs. / pallet	100000
Minimum order qty.	2000
Depth of case	458
Width of case	313
Height of case	180

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Micro tube 1.5ml SafeSeal

Page 2

Case volume	0,0258 cbm
Weight of case	2,94 KG
EAN of inner box	4038917118293
EAN of outer box	4038917078320

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2016-05-23

This document was prepared by EDP support and is valid without signature.



SARSTEDT