



## Product description

Order number	70.1189.105
Product description	Pipette tip, 1 - 200 µl, neutral, filling level rings, extra long, universally suitable for 2 - 200 µl pipettes by Eppendorf, Gilson, Finnpipette, Biohit, Rainin and Socorex and structurally identical designs, 96 pieces in the box

## Product characteristics

Colour code	Yellow
Graduation	Yes
Filter	No
Autoclavable (max. °C)	121 °C

## Size

Max. volume	200 µl
Min. volume	1 µl
Length of product	77 mm

## Material & colours

Product material	Polypropylene (PP)
------------------	--------------------

## Purity & certification

Product in accordance with MDD or IVDD	No marking necessary
Certificate of conformity	No
CE certified	No CE marking
Lot or serial number	always batched

## Packaging

Type of packaging	Box
Pcs. / inner box	96
Pcs. / box	960
Pcs. / pallet	120960
Minimum order qty.	960
Case volume	0,0103 cbm
Weight of case	1,94 KG
EAN of inner box	4038917314497

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





EAN of outer box

4038917314480

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2016-05-31

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**