



Product description

Order number	11.2452.001
Product description	V-Monovette® Urine, 10 ml, Vacuum system for closed urine transfer, 100x15 mm, round base, paper label, 50/inner box sterile

Product characteristics

Type of collection	Vacuum
Label/ Print	With paper label
Colour of print/label	Yellow
Graduation	No
Cap	Membrane cap
Type of cap	Cap
Light protection	No
Application	Urine sample tubes
Example of use	Urine collection, transport and storage
Base shape	Round base
Recommended centrifugation	400 x g
Centrifugation temperature	20
Centrifugation time	5 min
Connection	V-Monovette®

Size

Sample volume	10 ml
Diameter	15 mm
Length including cap	106 mm
Length excluding cap	100 mm

Material & colours

Product material	Polyethylene terephthalate (PET)
Colour of product	Transparent
Cap material	High Density Polyethylene (HD-PE)
Colour of cap	Yellow
Membrane material	Rubber

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product in accordance with MDD or IVDD	In-Vitro Product
CE certified	CE
Sterilisation	Gamma irradiation
Lot or serial number	always batched

Packaging

Type of packaging	50 pcs. / box
Pcs. / inner box	50
Pcs. / box	500
Pcs. / pallet	34000
Minimum order qty.	500
Depth of case	428
Width of case	398
Height of case	135
Case volume	0,023 cbm
Weight of case	5,19 KG
EAN of inner box	4038917244855
EAN of outer box	4038917244848

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

