








Sensor Accessories				
Masimo PN	Description	UM	Product Image	Notes
1596	<b>CleanShield® Multisite Wrap</b> For use with YI Sensors Adult/Pediatric/Neonatal 100/box, Non-sterile This product does not contain natural rubber latex.	EA		
1597	<b>Standard Multisite Wrap</b> For use with YI Sensors Adult/Pediatric/Neonatal 100/box, non-sterile This product does not contain natural rubber latex.	EA		
1598	<b>Standard Petite Wrap</b> For use with YI Sensors Adult slender digit/Ped/Neonatal adhesive 100/box, non-sterile This product does not contain natural rubber latex.	EA		
1600	<b>Accessory Adhesive Squares</b> For use with YI Sensors 144/pack, non-sterile	EA		
1903	<b>Transflectance Sensor Adhesive Pad</b> For use with TF-I Sensors 20/pack, non-sterile	EA		
16179	<b>Tele Plast</b> Philips SpO2 Connector Plug. Compatible with Philips MMS module.	EA		
1602	<b>Foam Wraps for LNOP®/LNCS®/M-LNCS™ YI Sensors</b> For use with YI Sensors 12/pack, non-sterile Washable This product does not contain natural rubber latex.	EA		
1601	<b>Accessory Clothing Clips</b> For securing reusable sensor cables to patient clothing 5/pack, non-sterile	EA	