

Enk Oxygen Flow Modulator Set

Intended for emergency transtracheal catheter ventilation when conventional ventilation by mask or endotracheal tube cannot be performed.



FEATURES and BENEFITS

- The set provides clinicians with a compact solution for quickly establishing an airway and using jet ventilation to stabilise a patient in an emergent situation.
- The wire-reinforced catheter is designed to resist kinking and prevent the airway device from collapsing after placement.
- The catheter connects to Luer lock adapters to allow oxygen delivery.

COMPONENTS

- Enk Oxygen Flow Modulator
- Emergency Transtracheal Airway Catheter
- Connecting tube
- Syringe

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter ID mm	Catheter Needle gage	Connecting Tube Length cm
G12398	C-EFMS-100	6	7.5	2	15	210

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.