



Product description

Order number	72.704.001
Product description	SafeSeal Micro tube 0.5 ml with attached lid with a wide hinge and retaining cams, PP, brown, 500 pieces packed in the bag

Product characteristics

Graduation	Yes
Cap	SafeSeal
Recommended centrifugation	30000 x g

Size

Sample volume	0,5 ml
Diameter	7,8 mm
Length of product	30 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	Brown
Cap material	Polypropylene (PP)

Purity & certification

Product in accordance with MDD or IVDD	No marking necessary
CE certified	No CE marking
Lot or serial number	not always batched

Packaging

Pcs. / inner box	500
Pcs. / box	2000
Pcs. / pallet	156000
Minimum order qty.	2000
Depth of case	388
Width of case	358
Height of case	135
Case volume	0,0188 cbm
Weight of case	1,61 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
SafeSeal tube 0.5ml, brown

Page 2

EAN of inner box 4038917227520
EAN of outer box 4038917227513

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2016-05-23

This document was prepared by EDP support and is valid without signature.



SARSTEDT