

# BP Cuff Adapters

## Connecting BP Cuff to DINACLICK and Neo-Snap Connectors



**Y, 1 Male Bayonet to DINACLICK\* Female Recept**

Compatibility: Philips/HP, Siemens, Datascope, Colin

Part # 2059291-001 (10/Pk)



**Y, 1 Female Submin to DINACLICK Female Recept**

Compatibility: Protocol

Part # 2059291-002 (10/Pk)



**Y, T Male Submin to DINACLICK Female Recept**

Compatibility: Nihon Kohden

Part # 2059291-003 (10/Pk)



**Y, 1 Unfitted to DINACLICK Female Recept**

Compatibility: Unfitted

Part # 2059291-004 (10/Pk)



**X, 2 Male Bayonet to DINACLICK Female Recept**

Compatibility: User Modified

Part # 2059292-001 (10/Pk)



**X, 2 Female Submin to DINACLICK Female Recept**

Compatibility: GE Marquette

Part # 2059292-002 (10/Pk)



**X, 2 Mated Submin to DINACLICK Female Recept**

Compatibility: GE Datex-Ohmeda

Part # 2059292-003 (10/Pk)



**X, 2 Male Screw to DINACLICK Female Recept**

Compatibility: DINAMAP

Part # 2059292-004 (10/Pk)



**Y, 1 Male Slip Luer to Neo-Snap\* Male Plug**

Compatibility: Philips/HP, Siemens, Datascope, Spacelabs

Part # 2059290-001 (10/Pk)



**Male Micro Plug to 2 Neo-Snap Male Plug**

Compatibility: Spacelabs

Part # 2059290-002 (10/Pk)



**Male Micro Plug to 1 Neo-Snap Male Plug**

Compatibility: Spacelabs

Part # 2059297-001 (10/Pk)