

## Specimen container, conical, with screw caps

Leverantör: VWR Collection



PP

General specimen containers. 120 ml capacity.

Complete range available with different cap colours

Sterile as well as non sterile, individually packed or bulk

With frosted labelling area

VWR art.nr	Antal	Typ	Lock, färg	Förp.	Version	Förp.
216-0861		Sterile	Röd	Bulk	Assembled cap	300
216-0856		Sterile	Röd	Individually	Assembled cap	250
216-0857		Sterile	Blå	Bulk	Assembled cap	300
216-0860		Sterile	Vit	Bulk	Assembled cap	300
216-0869		Sterile	Vit	Individually	Assembled cap	250
216-0862		Sterile	Grön	Bulk	Assembled cap	300
216-0855		Non sterile	Röd	50/påse	Separate cap	500
216-0859		Non sterile	Röd	Bulk	Separate cap	500
216-0864		Non sterile	Röd	Bulk	Assembled cap	300
216-0868		Non sterile	Röd	Individually	Assembled cap	250
216-0858		Non sterile	Blå	Bulk	Assembled cap	300
216-0863		Non sterile	Vit	Bulk	Assembled cap	300
216-0865		Non sterile	Grön	Bulk	Assembled cap	300
216-3949		Non sterile	Grön	Bulk	Separate cap	500
216-1185		Non sterile	—	Bulk	Without cap	1.000