

Specimen containers, conical, with screw caps

Leverantör: VWR Collection



PP.

- ▶ Complete range available with different cap colours
- ▶ Sterile as well as non sterile, individually packed or bulk
- ▶ Frosted labelling area

General specimen containers. 120 ml capacity.

ORDER

VWR art.nr	Antal	Typ	Lock, färg	Version	Förp.	Förp.
216-0861		Sterile	Röd	Assembled cap	Bulk	300
216-0856		Sterile	Röd	Assembled cap	Individually	250
216-0857		Sterile	Blå	Assembled cap	Bulk	300
216-0860		Sterile	Vit	Assembled cap	Bulk	300
216-0869		Sterile	Vit	Assembled cap	Individually	250
216-0862		Sterile	Grön	Assembled cap	Bulk	300
216-0855		Non sterile	Röd	Separate cap	50/påse	500
216-0859		Non sterile	Röd	Separate cap	Bulk	500
216-0864		Non sterile	Röd	Assembled cap	Bulk	300
216-0868		Non sterile	Röd	Assembled cap	Individually	250
216-0858		Non sterile	Blå	Assembled cap	Bulk	300
216-0863		Non sterile	Vit	Assembled cap	Bulk	300
216-0865		Non sterile	Grön	Assembled cap	Bulk	300
216-3949		Non sterile	Grön	Separate cap	Bulk	500
216-1185		Non sterile	—	Without cap	Bulk	1.000

Separate cap

VWR art.nr	Antal	Färg	För	Förp.
216-1288		Röd	Urine container, 120 ml	1.000