

## Coaxial Introducer Guide Wires



Size: 0.018" (0.46mm) x 60cm  
Core/Tip: Nitinol/Platinum Tip

**MAK001N-60**

### GUIDE WIRES (CONTAINED IN A HOOP AND INDIVIDUALLY POUCHED)

MAK001	-	-	0.018" (0.46 mm) x 40 cm	SS/SS	-	10
MAK001N	-	-	0.018" (0.46 mm) x 40 cm	NT/PT	-	10
MAK001N-60	-	-	0.018" (0.46 mm) x 60 cm	NT/PT	-	10
MAK001N-80	-	-	0.018" (0.46 mm) x 80 cm	NT/PT	-	10
MAKNTPD	-	-	0.018" (0.46 mm) x 40 cm	NT/PD	-	10
MAK005	-	-	0.018" (0.46 mm) x 40 cm	SS/PT	-	10
MAK005-60	-	-	0.018" (0.46 mm) x 60 cm	SS/PT	-	10
MAKSSPD	-	-	0.018" (0.46 mm) x 40 cm	SS/PD	-	10



PERIPHERAL INTERVENTION

# ACCESS