



Product description

Order number	20.1346.100
Product description	Microvette® 500 µl, Lithium Heparin Gel, green US code, print

Product characteristics

Type of collection	Capillary
Type of preparation	Gel/heparin
Colour code	US
Label/ Print	With printed label
Colour of print/label	Green
Cap	Quadruple-threaded cap
Application	Capillary blood collection
Base shape	Flat base
Recommended centrifugation	10000 x g
Centrifugation temperature	Room temperature
Centrifugation time	5 min
Shape of inner tube	Cylindrical

Size

Sample volume	500 µl
Diameter	10,8 mm
Length including cap	47,6 mm
Length excluding cap	46,6 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	Transparent
Cap material	High Density Polyethylene (HD-PE)
Colour of cap	Green

Purity & certification

Product in accordance with MDD or IVDD	In-Vitro Product
CE certified	CE
Lot or serial number	always batched

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Pcs. / inner box	100
Pcs. / box	1000
Pcs. / pallet	121000
Minimum order qty.	1000
Depth of case	338
Width of case	228
Height of case	160
Case volume	0,0123 cbm
Weight of case	2,67 kg
EAN of inner box	4038917127868
EAN of outer box	4038917028790

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

