



Product description

Order number	70.1186.210
Product description	Filter tip, 50 - 1.250 µl, neutral, filling level rings, extra long, Biosphere® certified, universally suitable for 1.000 µl pipettes and for 1.200 µl multi-channel pipettes by Eppendorf and Finnpipette, 96 pieces in the box

Product characteristics

Colour code	Blue
Graduation	Yes
Filter	Yes

Size

Max. volume	1250 µl
Min. volume	50 µl
Length of product	103 mm

Material & colours

Product material	Polypropylene (PP)
Filter material	Polyethylene (PE)

Purity & certification

Product in accordance with MDD or IVDD	No marking necessary
Certificate of conformity	No
CE certified	No CE marking
Purity standard	Biosphere®
Sterilisation	Gamma irradiation
Pyrogen-free	Pyrogen-free
Lot or serial number	always batched

Packaging

Type of packaging	Box
Pcs. / inner box	96
Pcs. / box	768
Pcs. / pallet	77568
Minimum order qty.	768

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Biosphere®Fil.Tip 1250µl, long

Page 2

Depth of box	150
Width of box	20
Height of box	60
Case volume	0,0152 cbm
Weight of case	2,08 kg
EAN of inner box	4038917126496
EAN of outer box	4038917080354

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2016-05-31

This document was prepared by EDP support and is valid without signature.



SARSTEDT