



### **Product description**

Order number	95.1006
Product description	Disposable tournistrip® tourniquet made of paper, latex-free, 200 pieces packed in the donation case

### **Product characteristics**

Design	Disposable tourniquet
Application	Venous blood collection
Example of use	Disposable tourniquet for the collection of venous blood
Suitable for	95.1005

### **Material & colours**

Product material	Special paper/plastic composite material
------------------	--

### **Purity & certification**

Product in accordance with MDD or IVDD	Medical Product
CE certified	Supplier indicated on the product
Lot or serial number	always batched

### **Packaging**

Type of packaging	On the roll
Pcs. / inner box	200
Pcs. / box	4800
Pcs. / pallet	201600
Minimum order qty.	200
Depth of case	388
Width of case	378
Height of case	235
Case volume	0,0345 cbm
Weight of case	12,43 KG
EAN of inner box	4038917351003
EAN of outer box	4038917350990

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

