

General Purpose - Transducer Covers

Latex Free Telescopically Folded Extended Length Covers

The covers are available of 91.5 cm and 147 cm.

These extended length covers offer additional sterile protection when performing routine diagnostic procedures as well as puncture and drainage procedures.

The telescopic-fold allows for easy application of gel, as well as easy application of the cover onto the transducer.

Sterile general purpose procedure kits include transducer cover, sterile gel packet and colored elastic bands.

Select covers offer a three dimensional "box end."

Flat-Folded Covers

This cover is a specially formulated latex-free material with scanning properties similar to human tissue. The soft, pliable material is comfortable for patients and especially popular for vascular access, nerve blocks and other sterile scanning procedures. Sterile general purpose procedure kits include transducer cover, gel packet and colored elastic bands. Select covers offer a three-dimensional "box end."

Features:

- Distortion-free imaging
- Latex-free
- Easy application
- Powder free



Flat-folded Covers

Ref. no.	Qty/box	Description
610-1269	24	Sterile high quality flexible cover 14 x 61 cm Includes transducer cover, sterile gel packet and colored elastic bands.

Telescopically Folded Extended Length

Ref. no.	Qty/box	Description
610-1270	24	Sterile high quality flexible cover 14 x 91.5 cm Includes transducer cover, sterile gel packet and colored elastic bands.
610-1271	24	Sterile high quality flexible cover 14 x 147 cm Includes transducer cover, sterile gel packet and colored elastic bands.
610-1272	24	Sterile high quality flexible cover 8.9 x 91.5 cm Includes transducer cover, sterile gel packet and colored elastic bands.
610-1273	24	Sterile high quality flexible cover 10.2 x 147 cm Includes transducer cover, sterile gel packet and colored elastic bands.



Sterile gel

All the boxes are included the AQUASONIC® 100 Ultrasound Transmissions gel. 20g sterilized foil pouches.

