

Instruments

Taps

311.320 Tap for Cortex Screws \varnothing 3.5 mm,
length 110/50 mm

Predrilling: \varnothing 2.5 mm



311.330 Tap for Cortex Screws \varnothing 3.5 mm,
length 180/110 mm

Predrilling: \varnothing 2.5 mm



311.340 Tap for Cancellous Bone Screws \varnothing 4.0 mm,
length 110/85 mm



2-flute Drill Bit for Quick Coupling

310.230 Drill Bit \varnothing 2.5 mm, length 180/155 mm,
2-flute, for Quick Coupling



310.250 Drill Bit \varnothing 2.5 mm, length 110/85 mm,
2-flute, for Quick Coupling



310.350 Drill Bit \varnothing 3.5 mm, length 110/85 mm,
2-flute, for Quick Coupling



310.370 Drill Bit \varnothing 3.5 mm, length 195/170 mm,
2-flute, for Quick Coupling



*for predrilling pilot holes for Schanz Screws
 \varnothing 5.0 and 6.0 mm*

2-flute Drill Bit for Jacobs Chuck

310.240 Drill Bit \varnothing 2.5 mm, length 95/80 mm,
2-flute, for Jacobs Chuck

