

TWEEZER® Bone Marrow Biopsy Needle

The TWEEZER® Bone Marrow Biopsy Needle has a specimen capturing device that is used to minimize the need for external manipulation of the needle and to help the user ensure recovery of adequate specimens.

The ergonomically designed “Twist-Lock” handle, combined with the stylet point, allows easy penetration of hard bone and keeps the handle assembly intact during bone penetration.

The TWEEZER® Bone Marrow Biopsy Needle includes the unique TWEEZER® Extraction Cannula, striving to produce a quality marrow sampling by entrapping the specimen in the needle without damaging it, thus eliminating the painful needle deflection technique.

Features:

- **Extra Sharp Stylet Point** designed for quick and easy penetration
- **Double Cutting Edges** on the tapered distal cannula for superior coring
- **Unique Extraction Cannula** captures and retains marrow specimens
- **Centimeter Markings** for estimated depth of penetration
- **Ergonomic Handle** designed to minimize hand and arm fatigue
- **Marked Probe** to verify specimen length before extraction
- **NeedleVISE™** protection included for securing sharps after use



TWEEZER® Bone Marrow Biopsy Needle

Item no.	Size	Qty/box
RTN-84-T	8G x 4 inch	10
RTN-86-T	8G x 6 inch	10
RTN-114-T	11G x 4 inch	10
RTN-116-T	11G x 6 inch	10



* Patent No. 5,522,398 & 5,634,473
NeedleVISE is a trademark of Atrion Medical Products.

DS-1-05 USA Version 8

TWEEZER® Bone Marrow Biopsy Trays

TWEEZER® Bone Marrow Biopsy Trays		
Item no.	Specification	Qty/box
BMET-8A-T	Includes 8G x 4" TWEEZER® Needle and 15 G x 2.8" Aspiration Needle	10
BMET-8-T	Includes 8G x 4" TWEEZER® Needle	10
BMET-11A-T	Includes 11G x 4" TWEEZER® Needle and 15 G x 2.8" Aspiration Needle	10
BMET-11-T	Includes 11G x 4" TWEEZER® Needle	10

TWEEZER® Bone Marrow Biopsy Tray Contents Include:

Qty	Description
1	20cc Luer Lock Syringe
1	5cc Luer Lock Syringe
1	19G x 1 ½" Filter Needle
1	21G x 1 ½" Needle
1	25G x 5/8" Needle
1	1% Lidocaine Hydrochloride Injection 5 ml
1	Lidocaine Insert
1	Mini Scalpel
3	10% Providone-Iodine USP Swabsticks
1	Towel
1	Luer Cap
5	3" x 3" Gauze Sponges
1	2" x 3" Bandage
1	Fenestrated Drape
1	24" x 24" CSR Wrap
1	Tweezer® Bone Marrow Biopsy Needle
1	Bone Marrow Aspiration Needle
1	Marked Obturator/Probe
2	Direction Inserts - Biopsy and Aspiration Needles
1	Probe Guide



DS-1-05 USA Version 8

