

# ETHICON

a Johnson & Johnson company



**VICRYL\***  
(Polyglactin 910) Suture

# VICRYL\*

(Polyglactin 910) Braided and Monofilament

## SUTURE

<b>Description:</b>	<b>Braided</b>	For easy handling and secure knot tying.
	<b>Coated</b>	For smooth passage through tissue and easy knot tie down.
	<b>Synthetic</b>	For minimal tissue reaction.
	<b>Absorbable</b>	By hydrolysis. Essentially complete in 56-70 days.

<b>Tensile strength:</b>	<b>Post implantation</b>	<b>Approximate % original strength remaining</b>
	14 days	75%
	21 days	50% (6-0 and larger) 40% (7-0 and smaller)
	28 days	25% (6-0 and larger)

**Color:** Violet or undyed.

**Range:** 10-0 to 2 (USP). Supplied as needled sutures and ligatures. Sizes 10-0 and 9-0 are monofilament in structure.

**Indications:** VICRYL\* sutures are intended for use in soft tissue approximation and/or ligation, including use in ophthalmic surgery, peripheral nerve adaptation and microsurgery for vessels less than 2 mm in diameter. Typical areas of use include:








- Fascia closure - due to the unique 28-day claim
- Subcutaneous fat
- Joint capsule
- Uterus

Full details in the Instruction for Use included in every package.







# VICRYL\* (Polyglactin 910) Braided or Monofilament Suture

## 1/2 Circle Taper Point Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 TF-6	6,5 mm	7-0	VIOLET	30		36	V1028H
 TF-1	10 mm	5-0	VIOLET	70		36	V133H
 RB-2	13 mm	6-0	VIOLET	70		36	PZ2571H
		6-0	UNDYED	45		12	W9981
		5-0	UNDYED	70	450 microns	36	V433H
 TF	13 mm	6-0	VIOLET	70		36	V926H
		5-0	VIOLET	70	350 microns	36	V925H
		4-0	VIOLET	70		36	V924H
 RB-1	17 mm	6-0	VIOLET	70		36	V302H
		5-0	VIOLET	70		36	V303H
		5-0	VIOLET	70	VISI-BLACK*	36	V3030H
		5-0	UNDYED	70		36	V213H
		4-0	VIOLET	70		36	V304H
		4-0	VIOLET	70	VISI-BLACK*	36	V3040H
		4-0	UNDYED	70		36	V214H
		4-0	VIOLET	8x45 C/R		12	V734D
		4-0	VIOLET	4x45 C/R		12	V735G
		4-0	VIOLET	8x70 C/R		12	V799D
		3-0	VIOLET	70		36	V305H
		3-0	UNDYED	70		36	V215H
		3-0	VIOLET	8x45 C/R		12	J713D
2-0	VIOLET	70		36	V306H		
 RB-1 RB-1	17 mm	3-0	VIOLET	70		36	V667H
 SH-1	22 mm	4-0	VIOLET	70		36	V310H
		4-0	VIOLET	70	VISI-BLACK*	36	V3100H
		4-0	UNDYED	70		36	V426H
		3-0	VIOLET	70		36	V311H
		3-0	VIOLET	70	VISI-BLACK*	36	V3110H
		3-0	UNDYED	70		36	V219H
		3-0	VIOLET	4x45 C/R		36	V782H
		2-0	VIOLET	70		36	V312H
		2-0	VIOLET	70	VISI-BLACK*	36	V3120H
		2-0	UNDYED	70		36	V220H
		2-0	VIOLET	8x45 C/R		24	V792E









# VICRYL\* (Polyglactin 910) Braided or Monofilament Suture

## 1/2 Circle Taper Point Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 SH	26 mm	4-0	VIOLET	70		36	V315H
		4-0	VIOLET	70	VISI-BLACK*	36	V3150H
		4-0	VIOLET	8x70 C/R		12	V783G
		3-0	VIOLET	70		36	V316H
		3-0	VIOLET	70	VISI-BLACK*	36	V3160H
		3-0	UNDYED	70		36	V416H
		3-0	UNDYED	70	VISI-BLACK*	36	V4416H
		3-0	VIOLET	90		36	V242H
		3-0	VIOLET	8x45 C/R		12	J774D
		3-0	VIOLET	8x45 C/R	VISI-BLACK*	24	V7740E
		3-0	VIOLET	8X70 C/R	VISI-BLACK*	24	V7840E
		2-0	VIOLET	70		36	V317H
		2-0	VIOLET	70	VISI-BLACK*	36	V3170H
		2-0	UNDYED	70		36	V224H
		2-0	VIOLET	90		36	V243H
		2-0	VIOLET	4X45 C/R		36	V779H
		2-0	VIOLET	8x45 C/R		12	J775D
		0	VIOLET	70		36	V318H
		0	UNDYED	70		36	V225H
 MH-1	31 mm	3-0	VIOLET	70		36	V245H
		2-0	VIOLET	70		36	V246H
		0	VIOLET	70		36	V247H
 MH	36 mm	4-0	VIOLET	90		36	PG2592H
		3-0	VIOLET	70		36	V322H
		3-0	VIOLET	90		36	V100H
		2-0	VIOLET	70		36	V323H
		2-0	UNDYED	70		36	V226H
		2-0	VIOLET	90		36	V101H
		0	VIOLET	70		36	V324H
		0	VIOLET	90		36	V102H
1	VIOLET	70		36	V325H		
 LH	40 mm	2-0	VIOLET	75		12	W9150






# VICRYL\* (Polyglactin 910) Braided or Monofilament Suture

## 1/2 Circle Taper Point Needles (Heavier wire diameter)

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 UCL	13 mm	3-0	VIOLET	70		36	V113H
		0	VIOLET	70		36	V115H
 UCLX	16 mm	3-0	VIOLET	70		36	V13H
		2-0	VIOLET	70		36	V14H
		0	VIOLET	70		36	V15H
		1	VIOLET	70		36	V16H
 CT-3	22 mm	3-0	VIOLET	70		36	V327H
		2-0	VIOLET	70		36	V328H
		0	VIOLET	70		36	V329H
 CT-2	26 mm	3-0	VIOLET	70		36	V332H
		2-0	VIOLET	70		36	V333H
		2-0	UNDYED	70		36	V269H
		2-0	VIOLET	4x45 C/R		24	V7642E
		0	VIOLET	70		36	V334H
		0	UNDYED	70		36	V270H
		0	VIOLET	8x45 C/R		24	V744E
		0	VIOLET	4x70 C/R		24	4PK2575YE
		1	VIOLET	70		36	V335H
2	VIOLET	70		36	V336H		
 MO-6	27 mm	1	VIOLET	90		36	V237H
 MO-5	31 mm	0	VIOLET	75		12	W9215
		1	VIOLET	75		12	W9216
 CT-1	36 mm	3-0	VIOLET	70		36	V338H
		3-0	UNDYED	90		36	V944H
		2-0	VIOLET	70		36	V339H
		2-0	VIOLET	90		36	V345H
		2-0	UNDYED	90		36	V945H
		2-0	UNDYED	8x70 C/R		12	JJ42D
		0	VIOLET	70		36	V340H
		0	VIOLET	90		36	V346H
		0	UNDYED	90		36	V946H
		0	VIOLET	8x45 C/R		12	J740D
		1	VIOLET	70		36	V341H
		1	VIOLET	90		36	V347H
		1	VIOLET	4x70 C/R		12	V736G
		2	VIOLET	90		36	V348H
2	UNDYED	90		36	V948H		
 MO-4	36 mm	0	VIOLET	70		36	V436H
		1	VIOLET	90		36	V255H





# VICRYL\* (Polyglactin 910) Braided or Monofilament Suture

## 1/2 Circle Taper Point Needles (Heavier wire diameter)


Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 CT	40 mm	3-0	VIOLET	70		36	V350H
		2-0	VIOLET	70		36	V351H
		2-0	VIOLET	8x45 C/R		12	8PM2544YG
		0	VIOLET	70		36	V352H
		0	VIOLET	90		36	V358H
		1	VIOLET	70		36	V353H
		1	VIOLET	90		36	V359H
		1	VIOLET	8x45 C/R		12	J753D
		2	VIOLET	90		36	V360H
	2	UNDYED	90		36	V960H	
 MO-45	45 mm	0	VIOLET	75		12	W9244
		1	VIOLET	75		12	W9245
		2	VIOLET	75		12	W9246
 CTX	48 mm	3-0	VIOLET	70		36	V362H
		2-0	VIOLET	70		36	V363H
		0	VIOLET	70		36	V364H
		0	VIOLET	90		36	V370H
		0	UNDYED	90		36	V170H
		0	VIOLET	8x45 C/R		12	J764D
		1	VIOLET	70		36	V365H
		1	VIOLET	90		36	V371H
		1	UNDYED	90		36	V366H
		1	VIOLET	8x45 C/R		12	J765D
		2	VIOLET	90		36	V372H
		2	UNDYED	90		36	V367H
		6	VIOLET	120		24	PN2512E
 TP-1	65 mm	2	VIOLET	90		12	V585G
 XLH	70 mm	2-0	VIOLET	70		12	V581G
		0	VIOLET	70		12	V582G

# VICRYL\* (Polyglactin 910) Braided or Monofilament Suture





## 1/2 Circle ETHIGUARD\* Blunt Point Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 MHD-1	31 mm	3-0	VIOLET	75		12	W9983
		2-0	VIOLET	75		12	W9984
 CTD-1	36 mm	2-0	VIOLET	75		12	W9985
		2-0	VIOLET	75	Heavy Needle	12	W9993
		0	VIOLET	75		12	W9986
		0	VIOLET	90	Heavy Needle	12	W9994
		1	VIOLET	90	Heavy Needle	12	W9995
 CTD	40 mm	0	VIOLET	75		12	W9988
		0	VIOLET	90	Heavy Needle	12	W9996
		1	VIOLET	90	Heavy Needle	12	W9997
 CTXD	48 mm	0	VIOLET	90		12	W9990
		1	VIOLET	90		12	W9991

## 3/8 Circle Blunt Point Needles



Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 BP-5	65 mm	1	VIOLET	100		12	W9391

## 5/8 Circle Taper Point Needles




Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 UR-6	26 mm	3-0	VIOLET	70		36	J601H
		2-0	VIOLET	70		36	V602H
		0	VIOLET	70		36	V603H
		1	VIOLET	70		36	VC2575H
 UR-6	26 mm	2-0	VIOLET	70		36	VCC2573H
 UR-5	36 mm	2-0	VIOLET	70		36	V375H
		0	VIOLET	70		36	V376H
		1	VIOLET	70		36	J377H
 UR-4	40 mm	2-0	VIOLET	70		36	V380H
		0	VIOLET	70		36	J381H

# VICRYL\* (Polyglactin 910) Braided or Monofilament Suture






## Fish Hook Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 PB-30	30 mm	2-0	VIOLET	70	Taper point	36	DF2570H
 V-38	31 mm	1	VIOLET	75	TAPERCUT*	12	W9335

## Straight Taper Point Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 ST-4	19 mm	3-0	VIOLET	70		36	BA2574H
 ST-4	19 mm	4-0	VIOLET	35		36	V620H
 ST	48 mm	2-0	VIOLET	90		36	V399H




## 1/2 Circle TAPERCUT\* Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 V-5	17 mm	4-0	VIOLET	70		12	V994G
		4-0	VIOLET	70		36	V994H
		4-0	UNDYED	70		36	V230H
		3-0	VIOLET	70		36	V276H
 V-6	22 mm	4-0	UNDYED	70		36	V231H
 V-7	26 mm	3-0	VIOLET	70		12	V998G
		3-0	VIOLET	70		36	V998H
		2-0	VIOLET	70		36	V999H
		2-0	VIOLET	6x45 C/R		24	6OB2545YE
		0	VIOLET	70		36	V279H
 V-30	31 mm	2-0	VIOLET	75		12	W9360
		0	VIOLET	75		12	W9361
		1	VIOLET	75		12	W9362
 V-34	36 mm	3-0	VIOLET	90		36	V516H
		2-0	VIOLET	90		36	V517H
		0	VIOLET	90		36	V518H
		1	VIOLET	90		36	V519H






# VICRYL\* (Polyglactin 910) Braided or Monofilament Suture



## 1/2 Circle TAPERCUT\* Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 V-37	40 MM	2-0	UNDYED	75	Slimmed	12	W9900
		2-0	UNDYED	75		12	W9904
		2-0	VIOLET	90		12	W9466
		0	VIOLET	90		36	VS4000H
 V-39	45 mm	0	VIOLET	75		12	W9376
		1	VIOLET	75		12	W9377
 V-40	48 mm	2-0	VIOLET	90		36	VS5020H
		1	VIOLET	150	Loop	24	V829E
		6	VIOLET	75		12	V140G


## 3/8 Circle TAPERCUT\* Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 CC-1	13 mm	6-0	UNDYED	45	MULTIPASS*	12	JV1727
 CC-1	13 mm	6-0	VIOLET	75		36	JV1012
		5-0	VIOLET	75		36	JV1013
 V-4	17 mm	5-0	UNDYED	75		12	W9824
		4-0	UNDYED	70		12	V240G
		3-0	UNDYED	75		12	W9826

## 5/8 Circle TAPERCUT\* Needles







Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 UV-10	10 mm	6-0	VIOLET	30		12	W9172
 UV-17	17 mm	4-0	UNDYED	75		12	W9820

## Compound Curve TAPERCUT\* Needles



Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 DA	25 mm	4-0	VIOLET	45	VISI-BLACK*	12	18531G

# VICRYL\* (Polyglactin 910) Braided or Monofilament Suture

## 1/2 Circle Reverse Cutting Needles









Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
	22 mm	3-0	VIOLET	70		36	V460H
		3-0	UNDYED	70		36	ME2673H
		2-0	VIOLET	70		36	V461H
		2-0	UNDYED	75		36	JV459
	22 mm	0	UNDYED	70		36	J694H
		2	VIOLET	75		12	W9287
	26 mm	3-0	VIOLET	70		36	MF2572H
		2-0	VIOLET	75		36	JV969
		2-0	UNDYED	70		36	V869H
		0	VIOLET	70		36	MF2574H
		0	VIOLET	4x45 C/R		24	V870E
		1	UNDYED	70		36	MF2675H
	36 mm	3-0	VIOLET	70		36	V465H
		2-0	VIOLET	70		36	V466H
		2-0	VIOLET	90		36	V472H
		2-0	UNDYED	90		36	JV1472
		0	VIOLET	70		36	V467H
		0	VIOLET	90		36	V473H
		0	VIOLET	4x45 C/R		12	V475G
		1	VIOLET	70		36	V468H
		1	VIOLET	90		36	V474H
		2	VIOLET	90		36	V482H
	40 mm	2-0	VIOLET	90		36	JV484
		0	VIOLET	70		36	V479H
		0	VIOLET	90		36	V485H
		1	VIOLET	70		36	V480H
		1	VIOLET	90		36	V486H
	48 mm	0	VIOLET	90		36	V1057H
		1	VIOLET	90		36	V1058H
		2	VIOLET	90		36	V1059H

## 1/2 Circle Reverse Cutting PRIME\* Needles

Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
	8 mm	5-0	UNDYED	45		36	V505H
	16 mm	5-0	UNDYED	45	Compound curve	12	W9436
		4-0	UNDYED	45		12	J507G
		4-0	UNDYED	45	Compound curve	12	W9437








# VICRYL\* (Polyglactin 910) Braided or Monofilament Suture

## 3/8 Circle Reverse Cutting PRIME\* Needles





Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 P-1	11 mm	6-0	UNDYED	45	MULTIPASS*	36	MPV489H
		6-0	UNDYED	45		36	V489H
		5-0	UNDYED	45	MULTIPASS*	36	MPV490H
		5-0	UNDYED	45		36	V490H
 P-3	13 mm	6-0	UNDYED	45	MULTIPASS*	36	JV492H
		6-0	UNDYED	45		36	V492H
		5-0	UNDYED	45	MULTIPASS*	36	MPV493H
		5-0	UNDYED	45		36	V493H
		4-0	UNDYED	45	MULTIPASS*	36	MPV494H
		4-0	UNDYED	45		36	V494H
 PS-3	16 mm	5-0	UNDYED	45	MULTIPASS*	36	MPV500H
		4-0	UNDYED	45		24	W9506T
 PS-2	19 mm	5-0	UNDYED	45	MULTIPASS*	36	MPV495H
		4-0	UNDYED	45	MULTIPASS*	36	MPV496H
		4-0	UNDYED	45		24	W9515T
		3-0	UNDYED	45	MULTIPASS*	36	MPV497H
		3-0	UNDYED	45		36	V497H
 PS-2 S	19 mm	4-0	UNDYED	45	Slimmed	36	V496H
 PS-1	24 mm	4-0	UNDYED	70		36	JE2672H
		3-0	UNDYED	45	MULTIPASS*	36	MPV683H
		3-0	UNDYED	45		36	V683H
		2-0	UNDYED	45	MULTIPASS*	36	MPV684H
		2-0	UNDYED	45		36	V684H
 PS	26 mm	3-0	UNDYED	70	MULTIPASS*	36	MPV9526H
		2-0	UNDYED	75		24	W9527T
 PSL	30 mm	2-0	VIOLET	70		36	V686H

# VICRYL\* (Polyglactin 910) Braided or Monofilament Suture

## 3/8 Circle Reverse Cutting Needles




Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
	13 mm	6-0	VIOLET	45		36	J384H
		5-0	VIOLET	45		36	V385H
		4-0	VIOLET	45		36	V386H
	13 mm	4-0	VIOLET	45		12	JV559
	16 mm	6-0	VIOLET	45		36	V387H
		5-0	VIOLET	45		36	V388H
		4-0	VIOLET	75	CONV. CUTTING	12	W9443
	19 mm	5-0	UNDYED	45		36	V291H
		5-0	VIOLET	45		36	V391H
		5-0	UNDYED	70		36	V421H
		4-0	UNDYED	45		36	V292ZH
		4-0	VIOLET	45		12	V392ZG
		3-0	UNDYED	45		36	V293H
		3-0	VIOLET	45		12	V393G
		3-0	VIOLET	45		36	V393H
		3-0	UNDYED	70		36	V423H
		3-0	VIOLET	70		36	V398H
	19 mm	4-0	UNDYED	45	Slimmed	36	V292H
		4-0	VIOLET	45	Slimmed	36	V392H
		4-0	UNDYED	70	Slimmed	36	V422H
		4-0	VIOLET	70	Slimmed	36	V397H
	24 mm	4-0	UNDYED	70		36	V441H
		4-0	VIOLET	70		36	V451H
		3-0	UNDYED	70		36	V442H
		3-0	VIOLET	70		36	V452H
		2-0	UNDYED	70		36	V443H
		2-0	VIOLET	70		36	V453H
	30 mm	3-0	UNDYED	70		36	V588H
		3-0	VIOLET	70		36	V585H
		2-0	UNDYED	75		36	JV589
		2-0	VIOLET	70		36	V586H
		0	UNDYED	70		36	V520H
		0	VIOLET	70		36	V587H

## 3/8 Circle Conventional Cutting PRIME\* Needles



Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
	16 mm	5-0	UNDYED	45		24	W9505T
		4-0	UNDYED	45	MULTIPASS*	36	MPV9506H
	19 mm	4-0	UNDYED	75		24	W9570T
		3-0	UNDYED	75		24	W9571T
	26 mm	4-0	UNDYED	70	MULTIPASS*	36	MPV9580H
		4-0	UNDYED	75		24	W9580T
		3-0	UNDYED	75		24	W9581T
	36 mm	3-0	UNDYED	75		24	W9531T

# VICRYL\* (Polyglactin 910) Braided or Monofilament Suture



## Straight Cutting Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
 SC-1 	13 mm	4-0	VIOLET	45	36	V383H
 KS	60 mm	4-0	UNDYED	70	36	J662H
		3-0	VIOLET	70	36	V523H
		3-0	UNDYED	70	36	V663H
		2-0	VIOLET	75	12	W9418
		2-0	UNDYED	75	12	W9718

## 1/2 Circle MICROPOINT\* Spatula

Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code
 GS-14 GS-14	5,5 mm	10-0	VIOLET	30	12	V450G
		9-0	VIOLET	30	12	V449G
		8-0	VIOLET	20	12	V574G
		7-0	VIOLET	30	36	ENN2534H
 TG175-8	7 mm	7-0	VIOLET	45	12	W9565

## 3/8 Circle MICROPOINT\* Spatula




Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 GS-9	6,5 mm	8-0	VIOLET	30	12	V120G	
		7-0	VIOLET	30	12	V119G	
 GS-9 GS-9	6,5 mm	10-0	VIOLET	30	Monofilament	12	V440G
		9-0	VIOLET	30		12	V549G
		9-0	VIOLET	30	Monofilament	12	V439G
		8-0	VIOLET	20		12	V547G
		8-0	VIOLET	30		12	V548G
		7-0	VIOLET	45		12	V546G

## 1/4 Circle MICROPOINT\* Spatula Needles

Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code
 TG100-8	6,5 mm	6-0	VIOLET	45	12	V544G

# VICRYL\* (Polyglactin 910) Braided or Monofilament Suture

## 3/8 Circle MICROPPOINT\* Reverse Cutting Needles

Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 TG140-4	6,2 mm	10-0	VIOLET	30	Monofilament	12	EMM2830G
 G-6	8 mm	6-0	VIOLET	45		12	V551G
 G-1	11 mm	6-0	VIOLET	45		12	W9537
 G-3	13 mm	5-0	VIOLET	45		12	V557G


## 3/8 Circle CS ULTIMA\* Spatula

Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 CS-140-6	6,2 mm	10-0	VIOLET	15		12	V429G


## 1/8 Circle CS ULTIMA\* Spatula

Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 CS44-6	6,5 mm	10-0	VIOLET	10	Monofilament	12	V965

## 3/8 Circle Spatula Needles




Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 SPAT-10	11 mm	5-0	UNDYED	45		12	W9761
		4-0	UNDYED	45		12	W9762

## 1/2 Circle Spatula Needles






Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 S-22	8 mm	6-0	VIOLET	45		12	W9566
		6-0	UNDYED	45		12	W9756
		5-0	VIOLET	45		12	W9567

# VICRYL\* (Polyglactin 910) Braided or Monofilament Suture

## 1/4 Circle Spatula Needles

Needle	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code	
 S-4	8 mm	4-0	VIOLET	45		12	V572G
 S-29	8 mm	6-0	VIOLET	45		12	J555G
 S-14	8 mm	6-0	VIOLET	45		12	V570G
		6-0	VIOLET	45	VISI-BLACK*	12	V575G
		6-0	UNDYED	45		12	V670G
		6-0	VIOLET	20		12	W9562
		5-0	VIOLET	45		12	V571G
		5-0	UNDYED	45		12	V671G

## 3/8 Circle Taper Point Needles

Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code	
 BV-6	4,8 mm	10-0	VIOLET	10	Monofilament	12	V970G
 BV-4	4,9 mm	10-0	VIOLET	15	Monofilament	12	V971G
 BV130-5	6,5 mm	9-0	VIOLET	15	Monofilament	12	V239G
		8-0	VIOLET	20		12	V400G
 BV-1	9,3 mm	8-0	VIOLET	30		12	W9577
		6-0	VIOLET	30		12	W9575
 BB	17 mm	4-0	VIOLET	75	VISI-BLACK*	36	JV2024
		3-0	VIOLET	75	VISI-BLACK*	36	JV2025

# VICRYL\* (Polyglactin 910) Braided or Monofilament Suture

## Ligatures

	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
LIGAPAK*	4-0	VIOLET	250	Reel	12	V1204G
	3-0	VIOLET	250	Reel	12	V1205G
	2-0	VIOLET	250	Reel	12	V1206G
	0	VIOLET	250	Reel	12	V1207G
	1	VIOLET	250	Reel	12	V1208G
	4-0	VIOLET	125	Reel	12	V204G
	3-0	VIOLET	125	Reel	12	V205G
	2-0	VIOLET	125	Reel	12	V206G
	0	VIOLET	125	Reel	12	V207G
SUTUPAK*	3-0	VIOLET	140		36	V614H
	2-0	VIOLET	140		36	V615H
	0	VIOLET	140		36	V616H
	1	VIOLET	140		36	V617H
	2	VIOLET	140		36	V618H
SUTUPAK*	4-0	VIOLET	2x70		24	V623E
	3-0	VIOLET	2x70		24	V624E
	2-0	VIOLET	2x70		24	V625E
	0	VIOLET	2x70		24	V626E
	1	VIOLET	2x70		36	V627H
	2	VIOLET	2x70		36	V628H
SUTUPAK*	5-0	VIOLET	3x45		36	V632H
	5-0	UNDYED	3x45		36	V642H
	4-0	VIOLET	3x45		36	V633H
	4-0	UNDYED	3x45		36	V643H
	3-0	VIOLET	3x45		36	V634H
	3-0	UNDYED	3x45		36	V644H
	2-0	VIOLET	3x45		36	V635H
	2-0	UNDYED	3x45		36	V645H
	0	VIOLET	3x45		24	V636E
	1	VIOLET	3x45		24	V637E
2	VIOLET	3x45		36	V638H	
SUTUPAK*	4-0	VIOLET	5x70		24	V1214E
	3-0	VIOLET	5x70		24	V1215E
	3-0	UNDYED	5x70		24	V963E
	2-0	VIOLET	5x70		24	V1216E
	2-0	UNDYED	5x70		25	V964E
	0	VIOLET	5x70		24	V1217E
	0	UNDYED	5x70		24	V965E
	1	VIOLET	5x70		36	V1218H
	2	VIOLET	5x70		36	V1219H
SUTUPAK*	5-0	VIOLET	6x45		36	V1223H
	4-0	VIOLET	6x45		36	V1224H
	3-0	VIOLET	6x45		36	V1225H
	3-0	UNDYED	6x45		36	V978H
	2-0	VIOLET	6x45		36	V1226H
	0	VIOLET	6x45		24	V1227E
	1	VIOLET	6x45		24	V1228E
	2	VIOLET	6x45		36	V1229H



# VICRYL\* (Polyglactin 910) Braided or Monofilament Suture

## Ligatures

	Gauge	Suture Color	Length (cm)	Other information	Number per box	Code
SUTUPAK*	4-0	VIOLET	12x45		36	V903H
	3-0	VIOLET	12x45		36	V904H
	3-0	UNDYED	12x45		36	V910H
	2-0	VIOLET	12x45		24	V905E
	2-0	UNDYED	12x45		24	V911E
	0	VIOLET	12x45		24	V906E
	0	UNDYED	12x45		24	V912E
	1	VIOLET	12x45		24	V907E