



Product description

Order number	14.1217
Product description	Haemo-Diff® , without smear edge, for producing blood smears, 500 piece(s)/case

Product characteristics

Design	without smear edge
Application	Blood smear

Size

Length of product	26.5 mm
-------------------	---------

Material & colours

Product material	Styrene-butadiene (SB)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
**Haemo-Diff® , without smear edge, for producing
blood smears**

Page 2

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 2000
Type of inner packaging case
Piece(s) / inner box 500
Piece(s) / case 2000
Piece(s) / pallet 648000
Depth of case 268 mm
Width of case 191 mm
Height of case 100 mm
Case volume 0.0051 cbm
Weight of product 0.0003 kg
Weight of case 0.81 kg
EAN of inner box 4038917085823
EAN case 4038917004459

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-11-04

This document was prepared by EDP support and is valid without signature.



SARSTEDT