



---

Personalized PDF Catalog  
Catalogue Generated den 4 januari 2020

---

# Pipettspetsar, epT.I.P.S.® Standard

Leverantör: EPPENDORF

These tips feature an ergonomically optimised cone geometry which provides an improved tip seal.



- ▶ Spetsar av hög kvalitet
- ▶ Universalspetsar
- ▶ Flexibel kona

epT.I.P.S., in combination with adequate pipettes, are intended for dispensing liquids in the volume range of 0,1 µL to 10 ml.










In vivo applications (applications in or on the human body) are not allowed. epT.I.P.S. and the appropriate pipettes only must be used by correspondingly trained specialists.

The user must have read thoroughly the operating instructions of the pipette and the instructions for use of the epT.I.P.S. and must have made himself familiar with the functioning of the pipette with the pipette tips.

Förpackning: Bulk i återförslutningsbara påsar.

## ORDER

VWR art.nr	Antal	Volym	Steril	Längd	Spets färg	Färgkod	Förp.	Förp
613-3500		0,1 - 10 µl	-	34 mm	Transparent	Mörkgrå	Bulk	1.00
613-3501		0,1 - 20 µl	-	40 mm	Transparent	Grå	Bulk	1.00
613-3502		0,5 - 20 µl	-	46 mm	Transparent	Ljusgrå	Bulk	1.00
613-3503		2 - 200 µl	-	53 mm	Yellow	Gul	Bulk	1.00
613-3111		2 - 200 µl	-	53 mm	Transparent	Gul	Bulk, 2 bags × 500 tips	1.00
613-3112		20 - 300 µl	-	55 mm	Transparent	Orange	Bulk, 2 bags × 500 tips	1.00

613-3504		20 - 300 µl	-	55 mm	Yellow	Orange	Bulk
613-3505		50 - 1000 µl	-	71 mm	Blue	Blå	Bulk
613-3113		50 - 1000 µl	-	71 mm	Transparent	Blå	Bulk, 2 bags × 500 tips
613-3506		50 - 1250 µl	-	76 mm	Transparent	Grön	Bulk
613-3106		50 - 1250 µl	-	103 mm	Transparent	Grön	Bulk, 2 bags × 500 tips
613-3508		100 - 5000 µl	-	120 mm	Transparent	Lila	Bulk
613-3507		250 - 2500 µl	-	115 mm	Transparent	Röd	Bulk
EPPE0030000.650		200 - 5000 µl	-	175 mm	Transparent	Lila	Bulk, 3 bags × 100 tips
613-4269		1000 - 10000 µl	-	165 mm	Transparent	Turkos	Bulk, L=165 mm
613-4270		1000 - 10000 µl	-	243 mm	Transparent	Turkos	Bulk, L=243 mm