

Svensk distributör:
 Entomed MedTech AB
 Storgatan 7
 233 31 Svedala
 040-48 11 40
 info@entomedmedtech.se



www.sibelmed.com

DISPOSABLE OTOSCOPE SPECULUM

MODEL	HG4	HG2	HM4	HM2	R4	R2	R4L3	R2L3
DIAMETER	4 mm	2.5 mm	4 mm	2.5 mm	4 mm	2.5 mm	4 mm	2.5 mm
COMPATIBLE	Heine (Diagnóstico, K100, Beta 100) Kawe (Eurolight)	Heine (Diagnóstico, K100, Beta 100) Kawe (Eurolight)	Heine (Kappa, Beta, Alpha, MiniLux, Mini FibrLux, Mini 2000, K180, Beta 200, Mini 3000) Kawe (piccolight) Riester (Ri-Mini, Pen-Scope) Welch Allyn (serie 524)	Heine (Kappa, Beta, Alpha, MiniLux, Mini FibrLux, Mini 2000, K180, Beta 200, Mini 3000) Kawe (piccolight) Riester (Ri-Mini, Pen-Scope) Welch Allyn (serie 524)	Riester (Ri-Scope L1 y L2, Ri-Star)	Riester (Ri-Scope L1 y L2, Ri-Star)	Riester (Ri-Scope L3)	Riester (Ri-Scope L3)

Available in 2 formats (from 1,000u.)



Dispenser Box of 250u.



20 Strips of 50u.

Use products certificated in accordance with 93/42/EEC Directive in order to guarantee the patient security as well as the quality of the test.