

Svensk distributör:  
 Entomed MedTech AB  
 Storgatan 7  
 233 31 Svedala  
 040-48 11 40  
 info@entomedmedtech.se



www.sibelmed.com

## DISPOSABLE OTOSCOPE SPECULUM

MODEL	HG4	HG2	HM4	HM2	R4	R2	R4L3	R2L3
<b>DIAMETER</b>	4 mm	2.5 mm	4 mm	2.5 mm	4 mm	2.5 mm	4 mm	2.5 mm
<b>COMPATIBLE</b>	<p><b>Heine</b> (Diagnóstico, K100, Beta 100)</p> <p><b>Kawe</b> (Eurolight)</p>	<p><b>Heine</b> (Diagnóstico, K100, Beta 100)</p> <p><b>Kawe</b> (Eurolight)</p>	<p><b>Heine</b> (Kappa, Beta, Alpha, MiniLux, Mini FibrLux, Mini 2000, K180, Beta 200, Mini 3000)</p> <p><b>Kawe</b> (piccolight)</p> <p><b>Riester</b> (Ri-Mini, Pen-Scope)</p> <p><b>Welch Allyn</b> (serie 524)</p>	<p><b>Heine</b> (Kappa, Beta, Alpha, MiniLux, Mini FibrLux, Mini 2000, K180, Beta 200, Mini 3000)</p> <p><b>Kawe</b> (piccolight)</p> <p><b>Riester</b> (Ri-Mini, Pen-Scope)</p> <p><b>Welch Allyn</b> (serie 524)</p>	<p><b>Riester</b> (Ri-Scope L1 y L2, Ri-Star)</p>	<p><b>Riester</b> (Ri-Scope L1 y L2, Ri-Star)</p>	<p><b>Riester</b> (Ri-Scope L3)</p>	<p><b>Riester</b> (Ri-Scope L3)</p>

Available in 2 formats (from 1,000u.)



Dispenser Box of 250u.



20 Strips of 50u.

**Use products certificated in accordance with 93/42/EEC Directive in order to guarantee the patient security as well as the quality of the test.**