


# 100 pcs. disposable sanitary headphone covers



Maintain sterile, clean and healthy conditions between those who publicly share the same headphones  
Enables headphone wearers to hear music while maintaining effective hygiene between multiple  
headphone and headset users.

Elastic - will stretch to cover outer portion of the shared headphones where they contact the head and ear.  
The large size can be used on ALL types of earmuff style headphones also MRI pneumatic headphones.

Headphone covers are bagged in quantities of 100 sterile stretch headphone covers in each sealed bag.  
Comfortable material is latex free and does not itch or cause any added pressure against the skin.

Eliminates need for using cleaning agents and disinfectant sprays on your headphones, keeping the  
headphones germ free without deteriorating the cushion material

Will stretch around the outer portion of earmuff style headphones (11 to 17 cm diameter earpiece).

**Item no.: S-2 100**

