

MFX2270



Total Intravenous Administration (TIVA) Extension Set

SKU/REF MFX2270

SKU/REF Name TIVA extension set 3 way

Description 1 Anti-Siphon Valve - 1 Back Check Valve - length 210 cm

Resources (1)

[📄 Infusion specialty disposables e-catalogue](#)

Specifications

Country of Manufacture	Malta
Classification	IIa
Quantity per Box	150 (50x3)
Tube Length (cm)	206
Total Length (cm)	210
Priming Volume (mL)	5.2
Tube Inner Diameter (mm)	1 (second tube 3)
Pump Use	Yes, maximum 2.5 bar
Shelf Life	5 Years
Sterilization Method	Ethylene Oxide
Polyvinyl chloride (PVC)	Yes
DEHP	DEHP-Free
Ways	3-Way
Back check valve	Yes
Anti-syphon valve	Yes
Clamp	Yes
Luer Lock	Rotating luer lock
Application Area	Anaesthesia