




# Stille Forceps



**STILLE**

# Stille Instruments -an extension of the surgeon's hand

Since 1841, Stille has been part of surgery of its finest from. From the introduction of the double action rongeur in the early 1870's, to breakthroughs in stainless steel (1923), Super Cut scissors (1983), and other innovations, Stille instruments are the choice of the world's finest surgeons.

## Tissue Forceps

### ADSON

- 102-9019- Very delicate tips, 1x2 fine teeth, 12 cm 4.75"
- 102-9022- Delicate tips, X-serrated, 1x2 fine teeth, 12 cm 4.75"
- 102-9022-15 Delicate tips, X-serrated, 1x2 fine teeth, 15 cm 6"



### ADSON-CALLISON

- 102-9021- Delicate tips, Tying platform, 1x2 fine teeth, 12 cm 4.75"
- 102-9021-1 Very delicate tips, Tying platform, 1x2 fine teeth, 12 cm 4.75"




### ADSON-BROWN

- 102-9023-12 Delicate tips, Side-grasping teeth, 12 cm 4.75"
- 102-9023-15 Delicate tips, Side-grasping teeth, 15 cm 6"



### GILLIES

- 102-9027- Delicate tips, Serrated, 1x2 fine teeth, 15.5 cm 6.25"



### STILLE

- 102-9028- Delicate tips, Serrated, 1x2 fine teeth, 15.5 cm 6.25"



### BISHOP-HARMON

- 102-9616- Very delicate tips, 1x2 fine teeth, 8.5 cm 3.4"



### ADLERQREUTZ

- 102-9248-15 Wide tips, 4x5 fine teeth, 15 cm 6"
- 102-9248-20 Wide tips, 4x5 fine teeth, 20 cm 8"
- 102-9248-25 Wide tips, 4x5 fine teeth, 25 cm 10"



### WAUGH-BROPHY

- 102-9252-20 Delicate X-serrated, 1x2 fine teeth, 20 cm 8"
- 102-9252-25 Delicate X-serrated, 1x2 fine teeth, 25 cm 10"
- 102-9252-30 Delicate X-serrated, 1x2 fine teeth, 30 cm 12"
- 102-9252-33 Delicate X-serrated, 1x2 fine teeth, 33 cm 13"



### RUSSIAN-STILLE

- 102-9270-15 Rounded wide tips, multiple fine teeth, 15 cm 6"
- 102-9270-20 Rounded wide tips, multiple fine teeth, 20 cm 8"



### STILLE

- 102-9050-11 1 x 2 fine teeth, 11 cm 4.25"
- 102-9050-13 1 x 2 fine teeth, 13 cm 5"
- 102-9060-11 1 x 2 teeth, 11 cm 4.25"
- 102-9060-15 1 x 2 teeth, 15 cm 6"
- 102-9060-17 1 x 2 teeth, 17 cm 6.75"
- 102-9060-20 1 x 2 teeth, 20 cm 8"
- 102-9062-15 2 x 3 teeth, 15 cm 6"
- 102-9062-17 2 x 3 teeth, 17 cm 6.75"
- 102-9062-20 2 x 3 teeth, 20 cm 8"
- 102-9064-15 3 x 4 teeth, 15 cm 6"
- 102-9064-17 3 x 4 teeth, 17 cm 6.75"
- 102-9064-20 3 x 4 teeth, 20 cm 8"



Superior quality and design features....



## Dressing Forceps


### ADSON

- 102-9017- Very delicate tips, Serrated, 12 cm 4.75"
- 102-9018- Very delicate tips, X-serrated, 12 cm 4.75"
- 102-9020- Delicate tips, Serrated, 12 cm 4.75"
- 102-9020-1 Delicate tips, Diamond jaws, 12 cm 4.75"
- 102-9020-15 Delicate tips, Serrated, 15 cm 6"
- 102-9024- Delicate tips, TC jaws, 12 cm 4.75"
- 102-9029-15 Delicate tips, TC jaws, 15 cm 6"



### McINDOE

- 102-9025- Delicate tips, Serrated, 15.5 cm 6.25"
- 102-9025-1 Delicate tips, Diamond jaws, 15.5 cm 6"



### CUSHING

- 102-9606-17 Rounded tips, Diamond jaws, 17 cm 6.75"
- 102-9607-17 Rounded tips, Diamond jaws, angled 45°, 17 cm 6.75"



### BISHOP-HARMON

- 102-9615- Very delicate tips, Serrated, 8.5 cm 3.4"



### CRAFOORD-STILLE

- 102-9077-20 Rounded tips, Serrated, 20 cm 8"
- 102-9078-20 Rounded tips, Diamond jaws, 20 cm 8"



### WAUGH-BROPHY

- 102-9254- Delicate tips, X-serrated, 20 cm 8"
- 102-9255-20 Delicate tips, Diamond jaws, 20 cm 8"
- 102-9255-25 Delicate tips, Diamond jaws, 25 cm 10"



### STILLE

- 102-9030-13 Square tips, Serrated, 13 cm 5"
- 102-9030-15 Square tips, Serrated, 15 cm 6"
- 102-9030-17 Square tips, Serrated, 17 cm 6.75"
- 102-9030-20 Square tips, Serrated, 20 cm 8"
- 102-9030-25 Square tips, Serrated, 25 cm 10"
- 102-9040-13 Rounded tips, Serrated, Rounded upper end, 13 cm 5"
- 102-9040-15 Rounded tips, Serrated, Rounded upper end, 15 cm 6"
- 102-9040-17 Rounded tips, Serrated, Pointed upper end, 17 cm 6.75"
- 102-9040-20 Rounded tips, Serrated, Pointed upper end, 20 cm 8"
- 102-9040-25 Rounded tips, Serrated, Pointed upper end, 25 cm 10"
- 102-9040-30 Rounded tips, Serrated, Pointed upper end 30 cm 12"



## Insulated Diathermy Forceps

### BJÖRK-STILLE INSULATED DIATHERMY FORCEPS

- 102-9845-15 Serrated tips, 15 cm 6"
- 102-9845-20 Serrated tips, 20 cm 8"
- 102-9845-25 Serrated tips, 25 cm 10"



## STILLE MICRO

117-B00-	Straight tips, 16 cm 6.25"
117-B05-	Extra delicate straight tips, 16 cm 6.25"
117-B15-	Extra delicate curved tips, 16 cm 6.25"
117-B110-15	Straight tips, Tying platform, 15 cm 6"
117-B110-18	Straight tips, Tying platform, 18 cm 7"
117-B110-21	Straight tips, Tying platform, 21 cm 8.25"
117-B111-15	Curved tips, Tying platform, 15 cm 6"
117-B111-18	Curved tips, Tying platform, 18 cm 7"
117-B111-21	Curved tips, Tying platform, 21 cm 8.25"
117-B120-15	Straight tips, Diamond jaws, Tying platform, 15 cm 6"
117-B120-18	Straight tips, Diamond jaws, Tying platform, 18 cm 7"
117-B120-21	Straight tips, Diamond jaws, Tying platform, 21 cm 8.25"
117-B121-15	Curved tips, Diamond jaws, Tying platform, 15 cm 6"
117-B121-18	Curved tips, Diamond jaws, Tying platform, 18 cm 7"
117-B121-21	Curved tips, Diamond jaws, Tying platform, 21 cm 8.25"
117-B130-18	Straight tips, Diamond jaws, Tying platform, Counter-balanced, 18 cm 7"
117-B130-21	Straight tips, Diamond jaws, Tying platform, Counter-balanced, 21 cm 8.25"
117-B131-18	Curved tips, Diamond jaws, Tying platform, Counter-balanced, 18 cm 7"
117-B131-21	Curved tips, Diamond jaws, Tying platform, Counter-balanced, 21 cm 8.25"
117-B140-18	Straight ring tips 1.0 x 0.5 mm, Diamond jaws, Counter-balanced, 18 cm 7"
117-B140-21	Straight ring tips 1.0 x 0.5 mm, Diamond jaws, Counter-balanced, 21 cm 8.25"
117-B145-18	Straight ring tips 0.8 x 0.4 mm, Diamond jaws, Counter-balanced, 18 cm 7"
117-B145-21	Straight ring tips 0.8 x 0.4 mm, Diamond jaws, Counter-balanced, 21 cm 8.25"



 Straight

 Curved

 Ring tips


## STILLE-DeBAKEY VASCULAR FORCEPS

- 102-9300- 2 mm DeBakey tips, 15 cm 6"
- 102-9302- 2 mm DeBakey tips, 20 cm 8"
- 102-9304- 2 mm DeBakey tips, 25 cm 10"
- 102-9312- 2.7 mm DeBakey tips, 15 cm 6"
- 102-9314- 2.7 mm DeBakey tips, 20 cm 8"
- 102-9316- 2.7 mm DeBakey tips, 25 cm 10"
- 102-9318- 2.7 mm DeBakey tips, 30 cm 12"
- 102-9320- 3.2 mm DeBakey tips, 20 cm 8"
- 102-9322- 3.2 mm DeBakey tips, 25 cm 10"
- 102-9308- 2 mm DeBakey tips, Angled 30°, 20 cm 8"
- 102-9330- 2 mm DeBakey tips, Craford handle, 20 cm 8"



## STILLE IRIS FORCEPS

- 102-9460- Delicate tips, Serrated, 10 cm 4"
- 102-9470- Delicate tips, 1 x 2 fine teeth, 10 cm 4"



## OLIVECRONA-STILLE DURA FORCEPS

- 102-9244-15 Square tips, 2 x 3 fine teeth, 15 cm 6"
- 102-9244-20 Square tips, 2 x 3 fine teeth, 20 cm 8"
- 102-9244-25 Square tips, 2 x 3 fine teeth, 25 cm 10"



## BJÖRK-STILLE KNOT HOLDING FORCEPS

- 102-9846-20 Serrated tips, 20 cm 8"



## STILLE-MOBERG SUTURE FORCEPS

- 102-9229- 15 cm 6"



## STILLE FIXATION FORCEPS

- 102-9611- Fine multiple teeth, 4 mm tip, 11 cm 4.3"



• Photo shows this instrument. Photos not in actual size.



Stille Surgical AB

Sundbybergsvägen 1A SE-171 73 Solna Sweden

Phone:+46 8 588 580 00 Fax:+46 8 588 580 05

Email:info@stille.se www.stille.se