


Ragnell Stille dissektionssax | Ragnell Stille Dissecting Scissors

Dissecting scissors with flat, rounded tips. Often used in plastic surgery. Is also called Killner by some manufacturers. The short models are used for undermining in plastic surgery and the long models are used for further dissection.



			Plastic	Orthopaedic	Popular Model
101-8378	Straight	12.5 cm 5"	●	●	
101-8380-12	Curved	12.5 cm 5"	●	●	●
101-8380-15	Curved	15 cm 6"	●	●	
101-8380-18	Curved	18 cm 7"	●	●	●
101-8380-20	Curved	20 cm 8"	●	●	●
101-8380-23	Curved	23 cm 9"	●	●	

