



Product description

Order number	60.610.001
Product description	Screw cap tube, 10 ml, 92 x 15.3 mm, round base, PP, with assembled cap, sterile, 100 pcs./bag

Product characteristics

Graduation	No
Cap	Assembled
Type of cap	Screw cap
Base shape	Round base
Centrifugation max (RCF)	3000 x g

Size

Sample volume	10 ml
Diameter	15,3 mm
Length of product	92 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polypropylene (PP)
Colour of product	Transparent
Cap material	High Density Polyethylene (HD-PE)
Colour of cap	Transparent
Tube material	Polypropylene (PP)

Purity & certification

Satisfies the requirement	IATA
Product in accordance with MDD or IVDD	No marking necessary
CE certified	No CE marking
Purity standard	Sterile
Sterilisation	sterile
Lot or serial number	batched

Packaging

Minimum order qty.	500
Type of packaging	Bag
Pcs. / inner box	100
Pcs. / box	500
Pcs. / pallet	22000
Depth of case	568 mm
Width of case	348 mm
Height of case	155 mm
Case volume	0,0306 cbm
Weight of case	3,14 KG
EAN of inner box	4038917114639
EAN of outer box	4038917014335

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

