



Product description

Order number	51.1534.500
Product description	Salivette® to determine the cortisol from the saliva, with blue cap and blue standard label, 100 pieces packed stacked in the box

Product characteristics

Label/ Print	With paper label
Cap	Push cap
Design	With synthetic swab
Application	Saliva examination to measure the cortisol
Recommended centrifugation	1000 x g
Centrifugation temperature	Room temperature
Centrifugation time	2 min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Diameter	16,8 mm
Length including cap	97 mm
Length excluding cap	86 mm

Material & colours

Product material	Polypropylene (PP)/ Low Density Polyethylene (LD-PE)
Cap material	Low Density Polyethylene (LD-PE)

Purity & certification

Product in accordance with MDD or IVDD	IVD
CE certified	CE
Lot or serial number	batched

Packaging

Minimum order qty.	500
Type of packaging	Box
Pcs. / inner box	100
Pcs. / box	500
Pcs. / pallet	36000
Depth of case	533 mm
Width of case	263 mm
Height of case	140 mm
Case volume	0,0196 cbm
Weight of case	4,09 KG
EAN of inner box	4038917159364
EAN of outer box	4038917159357

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

