

HAND AND FINGER SPLINTS



Frog Splint



- An aluminum splint with malleable arms that fold to hold the finger in the desired alignment.

Size	Padded	Unpadded
Small	100001P	100001
Medium	100002P	100002
Large	100003P	100003

Karstam Splint



- An Aluminum splint used to treat drop finger injuries.
- Maintains inter-phalangeal positions for injuries to distal and medial joints.
- Available in three sizes, padded and unpadded.

Size	Padded	Unpadded
Small	100290P	100290
Medium	100291P	100291
Large	100292P	100292