

PRODUCT DATA SHEET

Esmark Bandage

Designed with excellent compression and elastic properties to facilitate application of tourniquet to allow for essentially bloodless surgery

- Advanced washing technology results in minimum powder
- Does not contain Latex, Nano-Silver additives or Phthalates
- Lot traceable



Part Number	Description	Case Quantity
AR-35309	Bandage Esmark 7.5cm x 275cm	20/Cs
AR-35409	Bandage Esmark 10cm x 275cm	20/Cs
AR-35609	Bandage Esmark 15cm x 275cm	20/Cs
AR-35612	Bandage Esmark 15cm x 365cm	20/Cs



Delivering Excellence in Custom Solutions

For more information please contact us at:

ArcRoyal®
Virginia Rd, Kells,
Co Meath,
Ireland.

T : +353 46 928 0100
F : +353 46 928 0110
E : info@arcroyal.ie
W : www.arcroyal.ie