

RENASYS Dressing Kits

Foam Dressing Kits with Soft Port

Product description	Product code
<p>RENASYS-F</p> <p>Small</p> <ul style="list-style-type: none">• Foam block 10cm x 8cm x 3cm• RENASYS Transparent Film• RENASYS Soft Port 	66801668
<p>Medium</p> <ul style="list-style-type: none">• Foam block 20cm x 12.5cm x 3cm• 2 RENASYS Transparent Films• RENASYS Soft Port 	66801669
<p>Large</p> <ul style="list-style-type: none">• Foam block 25cm x 15cm x 3cm• 3 RENASYS Transparent Films• RENASYS Soft Port 	66801670
<p>X-Large</p> <ul style="list-style-type: none">• Foam block 48cm x 41cm x 1.5cm• 6 RENASYS Transparent Films• RENASYS Soft Port 	66801671