

Product Name	Reservoir bag, 2L with anti occlusion cage mount, 22F neck	
Reference numbers	2820000	
Product Type	Reservoir bag	
Manufacturer	Intersurgical	
Materials	PP Chloropren	
Packaging	Individually packaged	
Unit Packaging	10	
Box quantity	2560	
Pallet quantity		
Instructions for Use	Instructions for use with product	NO
	Instructions for use with each box	NO
Certification	CE0120	Class 2A
Applicable Standards	BS EN ISO 13485: Medical devices — Quality management systems — Requirements for regulatory purposes	
	ISO 9001: Quality Management Systems - Requirements	
	ISO 15223-1: Medical Devices - Symbols to be Used with Medical Device Labels, Labelling, & Information to be Supplied	
	BS EN ISO 18190: Anaesthetic and respiratory equipment. General requirements for airways and related equipment	
	ISO 5362: Anaesthetic reservoir bags	
	BS EN ISO 5356-1: Anaesthetic and respiratory equipment. Conical connectors. Cones and sockets	
Sterile	Not supplied sterile	
Shelf life	5 years	
Disposal Methods	As hospital guidelines	
Re-use	ⓧ	Single use only

Product characteristics

Technical information

Bag capacity (litres)	2
Leakage (ml/min)	3.3
Pressure head (mbar)	38



**Material
Abbreviations**

ABS	Acrylonitrile Butadiene Styrene
AL	Aluminium
HDPE	High Density Polyethylene
LDPE	Low Density Polyethylene
NRL	Natural Rubber Latex
PC	Polycarbonate
PE	Polyethylene
PP	Polypropylene
PVC	Poly Vinyl Chloride
SB	Styrene Butadiene
Si	Silicone
SS	Stainless Steel
TPE	Thermoplastic Elastomer

