


RENASYS ordering information

RENASYS TOUCH NPWT System

RENASYS TOUCH Device



Product description	Product code
RENASYS TOUCH <ul style="list-style-type: none"> • Pump • Power cord • RENASYS TOUCH Clinician User Manual • RENASYS TOUCH Home Healthcare User Manual 	66801280

RENASYS TOUCH Canisters

Product description	Product code
RENASYS TOUCH 300ml canister with solidifier <ul style="list-style-type: none"> • Sealed canister • Solidifier • Overflow guard with silicone tubing • Canister tubing 	66801273
RENASYS TOUCH 300ml canister without solidifier <ul style="list-style-type: none"> • Sealed canister • Overflow guard with silicone tubing • Canister tubing 	66801275
RENASYS TOUCH 800ml canister with solidifier <ul style="list-style-type: none"> • Sealed canister • Solidifier • Overflow guard with silicone tubing • Canister tubing 	66801274

RENASYS TOUCH

Accessories

Product description	Product code
RENASYS TOUCH Carry straps 	66801276
RENASYS TOUCH Carry bag 	66801277
RENASYS TOUCH IV pole/Bed clamp	66801278
RENASYS TOUCH Transit case	66801279
RENASYS TOUCH Odor Filter	66801284
RENASYS TOUCH O-Ring	66801283

RENASYS MAX NPWT System

RENASYS MAX Devices

Product description	Product code
RENASYS EZ MAX (Professional healthcare facility model) <ul style="list-style-type: none"> • Pump • Canister holder kit • Power cord • RENASYS EZ MAX Clinician User Manual 	66801309
RENASYS EZ Max (Home healthcare model) <ul style="list-style-type: none"> • Pump • Canister holder kit • Power cord • RENASYS EZ MAX Clinician User Manual • RENASYS EZ MAX Home Healthcare User Manual 	66801310

RENASYS EZ PLUS/MAX Canisters

Product description	Product code
RENASYS EZ 800ml S-canister 10 kits per case <ul style="list-style-type: none"> • Sealed canister • Overflow guard with silicone tubing • Solidifier • Canister tubing 	800ml with solidifier 66800912
	800ml without solidifier 66801066
RENASYS EZ 250ml S-canister with solidifier 10 kits per case <ul style="list-style-type: none"> • Sealed canister • Solidifier • Overflow guard with silicone tubing • Canister tubing 	250ml with solidifier 66800913

RENASYS EZ PLUS/MAX Accessories

Product description	Product code
RENASYS EZ Transit case	66800587
RENASYS EZ Canister holder	66800060

RENASYS GO NPWT System





RENASYS GO Device

Product description	Product code
RENASYS GO <ul style="list-style-type: none"> • Pump • Power cord • Quick reference guide • Strap • RENASYS GO Clinician User Manual • RENASYS GO Home Healthcare User Manual 	66801496






RENASYS GO Canisters

Product description	Product code
RENASYS GO Large 750ml canister 5 per case <ul style="list-style-type: none"> • Sealed canister • Solidifier • Canister tubing 	750ml with solidifier 66800916
RENASYS GO 300ml canister 5 per case <ul style="list-style-type: none"> • Sealed canister • Solidifier • Overflow guard with silicone tubing • Canister tubing 	300ml with solidifier 66800914
	300ml without solidifier 66800915





RENASYS GO Accessories

Product description	Product code
RENASYS GO Carrying case 	66800162
RENASYS GO Carry strap 	66800163
RENASYS GO O Rings 	66800603
RENASYS GO Odour Filters 	66800061
RENASYS GO Transit case	66800586

RENASYS Dressing Kits

Foam Dressing Kits with Soft Port	
Product description	Product code
RENASYS-F Small <ul style="list-style-type: none"> Foam block 10cm x 8cm x 3cm RENASYS Transparent Film RENASYS Soft Port 	66801668
Medium <ul style="list-style-type: none"> Foam block 20cm x 12.5cm x 3cm 2 RENASYS Transparent Films RENASYS Soft Port 	66801669
Large <ul style="list-style-type: none"> Foam block 25cm x 15cm x 3cm 3 RENASYS Transparent Films RENASYS Soft Port 	66801670
X-Large <ul style="list-style-type: none"> Foam block 48cm x 41cm x 1.5cm 6 RENASYS Transparent Films RENASYS Soft Port 	66801671
RENASYS -AB Abdominal Dressing kit with Soft Port <ul style="list-style-type: none"> 2 Foam blocks 40.5cm x 27.5cm x 2.5cm Organ protection layer 80cm x 66cm 6 RENASYS Transparent Films RENASYS Soft Port 	66801672

RENASYS Dressing Kits

Gauze Dressing Kits with Soft Port	
Product description	Product code
RENASYS-G Small <ul style="list-style-type: none"> AMD Gauze Dressing 15cm x 17cm RENASYS Transparent Film RENASYS Non-adherent gauze RENASYS Soft Port SECURA® No-Sting Barrier Film wipe Saline Paper ruler 	66800933
Medium <ul style="list-style-type: none"> 2 AMD Gauze Dressing 15cm x 17cm RENASYS Transparent Film RENASYS Non-adherent gauze RENASYS Soft Port SECURA No-Sting Barrier Film wipe Saline Paper ruler 	66800934
Large <ul style="list-style-type: none"> AMD Gauze Roll 11cm x 4m 2 RENASYS Transparent Films 2 RENASYS Non-adherent gauzes RENASYS Soft Port SECURA No-Sting Barrier Film wipe Saline Paper ruler 	66800935
X-Large <ul style="list-style-type: none"> 2 AMD Gauze Roll 11cm x 4m 3 RENASYS Transparent Film 4 RENASYS Non-adherent gauzes RENASYS Soft Port 2 SECURA No-Sting Barrier Film Wipes 2 Salines Paper ruler 	66800936
Sterile Gauze Dressing Kit with Soft Port <ul style="list-style-type: none"> AMD Gauze Roll 11cm x 4m AMD Gauze Dressing 15cm x 17cm 3 RENASYS Transparent Films RENASYS Soft Port 	66801665

RENASYS Drain Kits

Gauze Drain Dressing Kits	
Product description	Product code
RENASYS-G Gauze Dressing kit with drain 5 kits per case <ul style="list-style-type: none"> 1 Drain (available in 4 configurations listed) 1 RENASYS Quick-Click Adapter 1 Non-Adherent Gauze 7.5 x 20cm 1 Kerlix AMD Gauze Roll 11cm x 4m 1 Kerlix AMD Gauze Sponge 15cm x 17cm 1 RENASYS Transparent Film 4 Strips waterproof tape Ostomy paste 1 SECURA No-Sting Skin Prep barrier film wipe Saline Paper ruler 	10Fr Round Drain Gauze dressing kit 66801255 10Fr 10mm Flat Drain Gauze dressing kit 66801256 15Fr Channel Drain Gauze dressing kit 66801257 19Fr Round Drain Gauze dressing kit 66801258
RENASYS Drain Kits	
Product description	Product code
Drain Accessory kits 10 kits per case <ul style="list-style-type: none"> 1 Drain (available in 4 configurations listed) 1 RENASYS Quick-Click Adapter 2 Strips waterproof tape Ostomy paste 	10Fr Round Drain 66801251 10mm Flat Drain 66801252 15Fr Channel Drain 66801253 19Fr Round Drain 66801254

RENASYS Accessories

Kit Accessories	
Product description	Product code
NPWT Antimicrobial Large Gauze Roll 10 pouches per case <ul style="list-style-type: none"> 5 rolls of gauze individually packaged per pouch 	66800391
RENASYS Transparent Film 10 pouches per case <ul style="list-style-type: none"> 10 Transparent Films 20cm x 30cm per pouch 	66800394
RENASYS Transparent Film X-Large 10 pouches per case <ul style="list-style-type: none"> 5 Transparent Films 40cm x 60cm per pouch 	66801666
RENASYS Y-Connector 10 pouches per case <ul style="list-style-type: none"> 1 Y-Connector per pouch 	66800971
RENASYS Soft Port Stand-Alone 5 pouches per case <ul style="list-style-type: none"> 1 Soft Port per pouch 	66801673
RENASYS Adhesive Gel Patch 5 boxes per case <ul style="list-style-type: none"> 10 gel patches with double sided silicone adhesive hydrogel 10cm x 7cm per box 	66801082