

## Enteral Syringe Nutrifit ENFit LDT

Low dose tip syringes for enteral nutrition which can only be used with ENFit feeding tubes or with any device featuring a male ENFit connection.

The feeding tubes are fitted with an ENFit connection compliant with the ISO standard 80369-3. Consequently, accidental misconnections between devices intended for enteral nutrition and devices intended for other indications (venous access, spinal, epidural, vesical, respiratory...) are avoided and the right administration route is respected.

ENFit syringes are fitted with a purple-coloured plunger enabling instant identification.

3 volumes of syringes are available: 1 ml, 2.5 ml and 5 ml.

