



Product description

Order number	86.1253.001
Product description	Serological pipette 5 ml, with cotton plug, colour code: blue, polystyrene, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 pc./paper/plastic packaging

Product characteristics

Colour code	Blue
Colour of print/label	Blue
Graduation	Yes
Design	plugged, individual sterile

Size

Sample volume	5 ml
Diameter	8 mm
Length of product	340 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polystyrene (PS)
Colour of product	Transparent

Purity & certification

Product in accordance with MDD or IVDD	No marking necessary
CE certified	No CE marking
Sterilisation	sterile
Pyrogen-free	Pyrogen-free
Lot or serial number	batched

Packaging

Minimum order qty.	500
Pcs. / inner box	50
Pcs. / box	500
Pcs. / pallet	27000
Depth of case	410
Width of case	245
Height of case	295
Case volume	0,0305 cbm
Weight of case	4,54 KG
EAN of inner box	4038917127240
EAN of outer box	4038917021104

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

