



### Product description

Order number	86.1254.001
Product description	Serological pipette 10 ml, with cotton plug, colour code: orange, polystyrene, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 pc./paper/plastic packaging

### Product characteristics

Colour code	Orange
Colour of print/label	Orange
Graduation	Yes
Design	plugged, individual sterile

### Size

Sample volume	10 ml
Diameter	10 mm
Length of product	339 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material	Polystyrene (PS)
Colour of product	Transparent

### **Purity & certification**

Product in accordance with MDD or IVDD	No marking necessary
CE certified	No CE marking
Sterilisation	sterile
Pyrogen-free	Pyrogen-free
Lot or serial number	batched

### **Packaging**

Minimum order qty.	500
Pcs. / inner box	50
Pcs. / box	500
Pcs. / pallet	27000
Depth of case	410 mm
Width of case	245 mm
Height of case	325 mm
Case volume	0,0305 cbm
Weight of case	5,74 KG
EAN of inner box	4038917127295
EAN of outer box	4038917021128

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

