



Product description

Order number	72.690.001
Product description	Reaction tubes, 1.5 ml, PP, with attached lid, with moulded graduation and frosted writing space, 500 pcs./bag

Product characteristics

Graduation	Yes
Cap	Attached
Recommended centrifugation	20000 x g

Size

Sample volume	1,5 ml
Diameter	10,8 mm
Length of product	39 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	Transparent
Cap material	Polypropylene (PP)
Colour of cap	Transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product in accordance with MDD or IVDD	No marking necessary
CE certified	No CE marking
Lot or serial number	non batched

Packaging

Minimum order qty.	5000
Pcs. / inner box	500
Pcs. / box	5000
Pcs. / pallet	180000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0,0475 cbm
Weight of case	5,64 KG
EAN of inner box	4038917161657
EAN of outer box	4038917076005

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

