



### **Product description**

Order number	72.695.500
Product description	SafeSeal reaction tube, 2 ml, PP, with attached lid with retaining cams, with moulded graduation and frosted writing space, 250 pcs./bag

### **Product characteristics**

Graduation	Yes
Cap	SafeSeal
Recommended centrifugation	25000 x g

### **Size**

Sample volume	2 ml
Diameter	10,8 mm
Length of product	40 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	Transparent
Cap material	Polypropylene (PP)
Colour of cap	Neutral

### **Purity & certification**

Product in accordance with MDD or IVDD	No marking necessary
CE certified	No CE marking
Lot or serial number	non batched

### **Packaging**

Minimum order qty.	2000
Pcs. / inner box	1000
Pcs. / box	2000
Pcs. / pallet	100000
Depth of case	458 mm
Width of case	313 mm
Height of case	180 mm
Case volume	0,0258 cbm
Weight of case	3,09 KG
EAN of inner box	4038917266109
EAN of outer box	4038917266093

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

