



Product description

Order number	62.547.254
Product description	Screw cap tube, 50 ml, 114 x 28 mm, conical base, PP, with printed writing space and graduation, with assembled red cap, sterile and non-pyrogenic/endotoxin-free, 25 pcs./bag

Product characteristics

Label/ Print	With printed label
Colour of print/label	Blue/white
Graduation	Yes
Cap	Assembled
Type of cap	Screw cap
Base shape	Conical base
Centrifugation max (RCF)	15500 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Sample volume	50 ml
Diameter	28 mm
Length of product	114 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	Transparent
Cap material	High Density Polyethylene (HD-PE)
Colour of cap	Red
Tube material	Polypropylene (PP)

Purity & certification

Satisfies the requirement	ADR
Product in accordance with MDD or IVDD	IVD
CE certified	CE
Purity standard	Come Grow With Us
Sterilisation	sterile
Pyrogen-free	Pyrogen-free
Lot or serial number	batched

Packaging

Minimum order qty.	300
Type of packaging	Bag
Pcs. / inner box	25
Pcs. / box	300
Pcs. / pallet	8400
Depth of case	588 mm
Width of case	398 mm
Height of case	255 mm
Case volume	0,0597 cbm
Weight of case	4,52 KG
EAN of inner box	4038917111898
EAN of outer box	4038917013420

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

