# Product Specification Tube 15ml, 120x17mm, PP

Page 1



## **Product description**

Order number 62.554.502

Product description Tube 15 ml, 120 x 17 mm, conical base, PP, with print,

assembled red cap, 50 sterile and pyrogen-free pieces

packed in the bag

#### **Product characteristics**

Label/ Print With printed label

Colour of print/label White Graduation Yes

Cap Assembled
Type of cap Screw cap
Base shape Conical base
Centrifugation max (RCF) 15500 x g

Centrifugation temperature 20

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is <u>confidential</u> and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





## Product Specification Tube 15ml, 120x17mm, PP

#### **Size**

Sample volume 15 ml
Diameter 17 mm
Length of product 120 mm

#### **Material & colours**

Product material Polypropylene (PP)

Colour of product Transparent

Cap material High Density Polyethylene (HD-PE)

Colour of cap Red

### **Purity & certification**

Satisfies the requirement IATA Product in accordance with IVD

MDD or IVDD

CE certified CE
Purity standard Sterile
Sterilisation sterile

Pyrogen-free Pyrogen-free Lot or serial number batched

# **Packaging**

Minimum order qty. 500 Type of packaging Bag Pcs. / inner box 50 Pcs. / box 500 Pcs. / pallet 20000 Depth of case 598 mm Width of case 398 mm Height of case 165 mm Case volume 0,0393 cbm Weight of case 3,75 KG

EAN of inner box 4038917113182 EAN of outer box 4038917013932

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is <u>confidential</u> and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

