


New Optiflow™+ tracheostomy and mask interfaces

Introducing the OPT970 and OPT980

OPT970 tracheostomy interface

The new OPT970 tracheostomy interface improves upon the current OPT570 and OPT870 interfaces, featuring a new sputum guard designed with stability and ease of use in mind.



OPT980 mask interface adapter

The new OPT980 mask interface adapter improves upon the current RT013 interface, for connecting to face and tracheostomy mask interfaces with a standard 22 mm medical taper.



Both new interfaces feature:

An improved Evaqua™ membrane

Only Fisher & Paykel Healthcare products have Evaqua technology, which reduces the formation of mobile condensate near the patient. Our continuous improvement has led to the newest generation of Evaqua being **stronger** and **more durable**, for clinician confidence when delivering the therapy; and **quieter**, for patient comfort when sleeping. Furthermore, the new Evaqua tubing is **50 mm longer**, providing patients with more freedom of movement.

A new connector

We have made changes to the connector that fits to the heated breathing tube, making it smaller and less obtrusive, while still ensuring that it provides a stable fit for delivering optimal therapy.

Please contact your local representative for further information about this or any other Fisher & Paykel Healthcare product or therapy.