

StatLock® IV Stabilisation Devices

- Designed to limit catheter pistoning and micro-motion associated with tape securement
- Proven to significantly reduce overall IV complications; including phlebitis, infiltration, occlusion and dislodgement^{1,7,9}
- Proven to significantly reduce unscheduled catheter restarts and lengthen IV dwell time^{1,7,8,9}
- Designed to increase patient comfort
- Reduces needlestick risks to healthcare workers³

An Evidence-Based Stabilisation Device That Meets The 2011 INS Standards.¹

Clinical Studies Show That StatLock® IV Stabilisation Devices

- Have an average dwell time of 98 hours⁸
- May save nurses 29 minutes per IV treatment¹⁰
- Can provide a hospital an estimated annual cost savings in excess of \$200,000⁷
- Reduce accidental needlesticks³



IV0525

Select Device Retainers

- “Over-the-top” retainer secures to the luer of the pre-packaged tubing extension set
- Compatible with most peripheral IV catheters



IV0502

Ultra Device Retainers

- “Over-the-top” retainer secures to the push tab of B. Braun® and BD® safety catheters



IV0501

IV Device Retainers

- “Lock-tight” retainer secures the catheter in place
- Accommodates a wide variety of peripheral IV catheters

Neonate Device Retainers


- Skin safe hydrocolloid adhesive combined with a Tyvek® spine for added strength and flexibility



VNN2

Huber Plus® Device Retainers

- Adhesive well pad and foam shims secure Huber Plus® needles



VHFP1

Device Anchor Pads

- Breathable, latex-free Tricot fabric provides added catheter security and patient comfort
- Closed cell foam designed to be gentle on the skin

PIV Complication ⁷	Tape	Statlock® Device	% Reduction	P Value
Patient Complications	48%	16%	↓67%	< 0.001
Unscheduled Restarts	71%	17%	↓76%	= 0.0001
Phlebitis	3.6%	0.7%	↓80%	< 0.001

Description	Product Code	Box Quantity
StatLock IV Select, Tricot Select Pad, Luer Retainer, 15cm Luer Lock Micro Ext. Set	IV0521	25
StatLock IV Select, Tricot Select Pad, Luer Retainer, 15cm Luer Lock Macro Ext. Set	IV0525	25
StatLock IV Select Paediatric, Tricot Select Pad, Luer Retainer, 10cm Luer Lock Micro Ext. Set	IV0502	25
StatLock IV Select Paediatric, Tricot Select Pad, Luer Retainer, 15cm Luer Lock Macro Ext. Set	IV0503	25
StatLock IV Ultra, Tricot Ultra Pad, Catheter Retainer	IV0520	50
StatLock IV Ultra Paediatric, Tricot Ultra Pad, Catheter Retainer	IV0500	50
StatLock IV Ultra Neonate, Ultra Pad, Blue Catheter Retainer	IV0501	50
StatLock IV, Adhesive Well Pad for Winged Catheters	IV0621	50
StatLock IV, Foam Pad, 30cm Luer Lock Micro Ext. Set	VIDML12	25
StatLock IV, Foam Pad, 15cm Luer Lock Micro Ext. Set	VIDML6	25
StatLock IV, Foam Pad, No Ext. Set	VIDMNS	50
StatLock Neonatal 2, Anchor Pads	VNN2	50
StatLock Huber Plus, Adhesive Well Pad, Foam Shims	VHFP1	50