



CLEAR  
SURGICAL  
SMOKE

**WITH  
CONFIDENCE.**



Valleylab™ Smoke Evacuation Pencil Portfolio








# Valleylab™ Smoke Evacuation Rocker Switch Pencil Edge™ Blade Electrode

	Medtronic Code	INS/ Old Code	Description	Length	Quantity
	SEP5000	INS-300-008 INS-500 INS-500-10	Smoke Evacuation Rocker Switch Pencil Edge™ Blade Electrode	10' (3 m)	20/case
	SEP5015	INS-515	Smoke Evacuation Rocker Switch Pencil Edge™ Blade Electrode	15' (4.6 m)	20/case
	SEP5000NSB	n/a	Smoke Evacuation Rocker Switch Pencil, NSB Edge™ Blade Electrode	10' (3 m)	30/case
	SEP5015NSB	n/a	Smoke Evacuation Rocker Switch Pencil, NSB Edge™ Blade Electrode	15' (4.6 m)	30/case
	SEP6000	INS-600	Telescoping Smoke Evacuation Rocker Switch Pencil Edge™ Blade Electrode	10' (3 m)	20/case
	SEP6015	INS-615	Telescoping Smoke Evacuation Rocker Switch Pencil Edge™ Blade Electrode	15' (4.6 m)	20/case
	SEP6000NSB	n/a	Telescoping Smoke Evacuation Rocker Switch Pencil, NSB Edge™ Blade Electrode	10' (3 m)	30/case
	SEP6015NSB	n/a	Telescoping Smoke Evacuation Rocker Switch Pencil, NSB Edge™ Blade Electrode	15' (4.6 m)	30/case













## Valleylab™ Smoke Evacuation Pencil Extension Tubes

	Medtronic Code	INS/ Old Code	Description	Length	Quantity
	SEA50595	INS-100-015	Smoke Evacuation Extension Tube For use with long electrodes	4" (10.16 cm)	80/case
	SEA54506	INS-410-080	Smoke Evacuation Extension Tube with Edge™ Blade Electrode For use with long electrodes	Tube 4" (10.16 cm) Electrode 6.5" (16.51 cm)	80/case





## Valleylab™ Tubing

	Medtronic Code	INS/ Old Code	Description	Length	Quantity
	SEA3705	INS-0316LEUR	Smoke Tubing	1/4" x 10" (0.6 cm x 3 m)	10/case
	SEA3710	INS-4400v	Smoke Tubing	3/8" x 10" (1 cm x 3 m)	10/case
	SEA3715	INS-102 INS-102S	Smoke Tubing	7/8" x 10" (2.2 cm x 3 m)	12/case
	SEA3720	INS-0316LEUR	Smoke Tubing	1/4" x 10" (0.6 cm x 3 m)	12/case
	SEA3725	n/a	Smoke Wand	7/8" x 8" (2.2 cm x 20.3 cm)	5/case


# Electrodes

	Medtronic Code	INS/ Old Code	Description	Length	Quantity
	E14504	n/a	Edge™ Blade Electrode	4" (10.16 cm)	50/case
	E14506	n/a	Edge™ Blade Electrode	6.5" (16.51 cm)	50/case
	E14526	n/a	Edge™ Needle Electrode	6.5" (16.51 cm)	50/case
	E14554	n/a	Edge™ Insulated Needle Electrode	4" (10.16 cm)	50/case
	E14556	n/a	Edge™ Insulated Blade Electrode	6.5" (16.51 cm)	50/case
	E1465	n/a	Edge™ Insulated Needle Electrode	2.8" (7.2 cm)	50/case
	E14654	n/a	Edge™ Insulated Needle Electrode	4" (10.16 cm)	50/case
	E14656	n/a	Edge™ Insulated Needle Electrode	6.5" (16.51 cm)	50/case
	E1502	n/a	Straight Extension	13 cm	1 each
	E1504	n/a	Straight Extension	34.3 mm	1 each
	E15516	n/a	Valleylab™ Blade Electrode	6.5" (16.51 cm)	50/case
	E15526	n/a	Valleylab™ Needle Electrode	6.5" (16.51 cm)	50/case

# Electrodes

	Medtronic Code	INS/ Old Code	Description	Length	Quantity
	E501564†	n/a	Ball Electrode 5 mm, for Standard SEP	2.5" (6.35 cm)	50/case
	E5015646†	n/a	Ball Electrode 5 mm, for Standard SEP	6.5" (16.51 cm)	50/case
	E601559†	n/a	Tungsten Loop Electrode 10 mm x 10 mm, for Telescoping	5" (13 cm)	10/case
	E601560†	n/a	Tungsten Loop Electrode 10 mm x 15 mm, for Telescoping	5" (13 cm)	10/case
	E601561†	n/a	Tungsten Loop Electrode 20 mm x 12 mm, for Telescoping	5" (13 cm)	10/case
	E601562†	n/a	Tungsten Loop Electrode 20 mm x 15 mm, for Telescoping	5" (13 cm)	10/case
	E601564†	n/a	Ball Electrode 5 mm, for Telescoping Smoke	4" (10.16 cm)	50/case
	E601651†	n/a	Straight Micro-needle Electrode for Telescoping Smoke	4" (10.16 cm)	50/case

## Compatible Electrodes for Valleylab™ Smoke Evacuation Pencil





†New



††Not compatible with an electrode less than 4 inches long and the SEA50595 extension tube.

†††Not compatible with an electrode less than 3.5 inches long.


## Valleylab™ Smoke Attachments

	Medtronic Code	INS/ Old Code	Description	Length	Quantity
	E3810	INS-3838 E3590	Valleylab™ Smoke Attachment	3/8 (1 cm) tubing, 10' (3 m)	20/case
	E3815	INS-3839	Valleylab™ Smoke Attachment	3/8 (1 cm) tubing, 15' (4.6 m)	20/case
	E3595	n/a	AccuVac™ Smoke Attachment Extender	n/a	10/case


## Valleylab™ Pencils with Smoke Attachments

	Medtronic Code	INS/ Old Code	Description	Length	Quantity
	E2350HS	n/a	Edge™ Rocker Pencil Holster Smoke Attachment	10' (3 m)	20/case
	E2350HSNSB	n/a	Edge™ Rocker Pencil Holster Smoke Attachment, NSB	10' (3 m)	50/case
	E2515HS	n/a	Rocker Pencil Holster Smoke Attachment	10' (3 m)	20/case
	E2515HSNSB	n/a	Rocker Pencil Holster Smoke Attachment, NSB	10' (3 m)	50/case
	E2450HS	n/a	Edge™ Pencil Holster Smoke Attachment	10' (3 m)	20/case
	E2450HSNSB	n/a	Edge™ Pencil Holster Smoke Attachment, NSB	10' (3 m)	50/case
	E2516HS	n/a	Button Pencil Holster Smoke Attachment	10' (3 m)	20/case





## Valleylab™ Fluid Trap

	Medtronic Code	INS/ Old Code	Description	Quantity
	RVFT10	INS-200 INS-200/EU INS-250	Fluid Trap	1 each

## Valleylab™ Filter

	Medtronic Code	INS/ Old Code	Description	Quantity
	SEA3700	n/a	ULPA Filter	2/case

## Valleylab™ Tubing Adapters

	Medtronic Code	INS/ Old Code	Description	Length	Quantity
	E3655	INS-011	OptiMumm™ Smoke Prefilter Reducer Adapter For use with non-Medtronic smoke evacuators.	7/8" to 3/8" (2.2 cm to 1 cm).	25/case
	E3656	n/a	Smoke Evacuation Pencil Adapter, IC Medical	3/8" (10 mm) to IC Medical Threaded Vacuum Connection	100/case
	E3660	INS-125	OptiMumm™ Smoke Prefilter Reducer Adapter For use with the OptiMumm™ smoke evacuator.	1 1/3" to 3/8" (3.4 cm to 1 cm).	25/case
	E3665	n/a	OptiMumm™ Smoke Prefilter Reducer Adapter For use with the OptiMumm™ smoke evacuator.	1 1/3" to 7/8" (3.4 cm to 2.2 cm).	25/case

Size	RapidVac™ Smoke Evacuator	RapidVac™ Fluid Trap	Opti- Mumm™ Smoke Adapter	Smoke Evac 200	IC Medical	Stryker Neptune™*	ConMed ClearView™*	ConMed Aer Defense™*	Buffalo ViroVac™*
E3655 (7/8" to 3/8")	No	No	No	Yes	No	No	No	No	No
E3660 (1 1/3" to 3/8")	No	No	Yes	No	No	No	No	No	No
E3665 (1 1/3" to 7/8")	No	No	Yes	No	No	No	No	No	No
SEP60XX SEP50XX E3810 E3815 with adapter	Yes	Yes	No	No	No	No	Yes	No	Yes
SEP60XX SEP50XX E3810 E3815 without adapter	No	Yes	No	No	No	Yes	Yes	No	Yes
E3656 (3/8" (10 mm) for IC Medical	No	No	No	No	Yes	No	No	No	No

For more information,  
call 800-772-8772  
or visit us at  
[medtronic.com/covidien](https://www.medtronic.com/covidien)