

# ONEMED

## PRODUKTBLAD

COVIDIEN - GIA6038S

### Wound Stapler GIA™ 60 Squeeze Handle Titanium Staples 3.8 mm



#### Features

- 3.8 mm / 60 mm
- White
- 60 Staples

#### Product Specifications

McKesson #	503002
Manufacturer #	GIA6038S
Brand	GIA™ 60
Manufacturer	Covidien
Application	Wound Stapler
Handle Type	Squeeze Handle
Staple Material	Titanium Staples
Staples Used	3.8 mm Staples
UNSPSC Code	42312009
Usage	Single Patient Use

#### Features

- 3.8 mm / 60 mm
- White
- 60 Staples
- The GIA™ 60 reloadable staplers place two double staggered rows of titanium staples and simultaneously cut and divide tissue between the two double rows
- Tear drop anvil / strongest beam available
- Lock lever on cartridge side
- Ergonomic inlay grips
- Rear hinge
- Adjustable firing knob
- Pinless for atraumatic placement
- Directional staple technology
- The GIA™ reloadable staplers should only be reloaded seven times per instrument for a total of eight firings
- WARNING: Cancer and Reproductive Harm - [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

Antal/avdelningsförpackning: 3  
Enhet ST